IKA-ULTRA-TURRAX® T 8

S 8 N - 5 N
S 8 N - 8 N

BETRIEBSANLEITUNG	DE	3
OPERATING INSTRUCTIONS	EN	10
MODE D'EMPLOI	FR	17

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DE

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
Contents

	Page
CE-Declaration of conformity	2
Safety instructions	10
Proper use	12
Unpacking	12
Useful information	12
Setting-up the instrument	12
Starting up	13
Drive unit ULTRA - TURRAX T8	13
Dispersing elements	13
Maintenance, Cleaning and Sterilisation	14
Accessories	15
Warranty	15
Technical data	15
Spare parts diagram/ - list of spare parts	16

Safety instructions

• Read the operating instructions in full before starting up and follow the safety instructions.

- In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.
- Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.
- Never operate the equipment with the dispersing tools rotating freely. Ensure that parts of the body, hair or items of clothing cannot be trapped by the rotating parts. Before switching on the dispersing instrument make sure that its shaft is immersed min. 10 mm in the medium to prevent the medium from splashing out.

- Do not use the appliance without a dispersing element.
- Vibrations can cause the turning handle on the dispersing tool to loosen. Check that the turning handle is secure and tighten if necessary.
- It is not permitted to use the dispersing instrument with a freely rotating dispersing element without a stirring vessel as any contact can cause serious injury. It is important to ensure that the dispersing element is always inside the stirring vessel when the appliance is switched on.
- The equipment is also suitable for manual operation.
- Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory.
- Please note that the dispersing element and the journal bearings can become extremely hot during use.
- When working with the dispersing unit, the user must select and wear his personal protective equipment according to the mixing hazard category. Defective or inappropriate protective equipment can expose the user to the risk of splashing or getting caught in the dispersing tool.
- Protect the dispersing element from any impact.
- Even small, undetectable damage creates unbalance and blocks the dispersing instrument. It also causes the dispersing element to overheat and inevitably overloads the drive.
- In the event of unbalance or unusual noises, switch off the appliance immediately in order to prevent unwanted resonance both in the appliance and the whole assembly. Any unbalance can damage the assembly, the appliance or destroy the stirring vessel, which, in turn, can seriously injure the user, as can the rotating dispersing element.
- Unusual noises indicate that the appliance and the dispersion tool is no longer working properly.
- Replace the dispersing element if you notice any unbalance or unusual noises. Stop using the appliance if there is no difference after you have changed the dispersing element and return it to the dealer or the manufacturer along with a description of the fault..
- Only dispersing elements approved by IKA may be used!

- Even in the case of repair work, the unit must only be opened by a qualified engineer. The mains plug must be pulled before opening the unit. The voltage-carrying components inside the unit may continue to carry a voltage for some time after the mains plug is pulled.
- With the hand enterprise and with the enterprise at the stand, it is to be made certain that safe conditions, a suitable anti-twist plate of the agitating container and a minimum distance of the dispersion tool are ensured to the agitating container soil.
- Do not cover the ventilation slots on the drive in order to ensure adequate cooling of the drive.
- The agitated vessels used for stirring have to be secured. Consider on a good stability of the entire structure.
- Furthermore, the dispersing unit may only be used to stir those materials or material mixtures that the user knows will not react dangerously to the extra energy produced by the dispersing. This also applies to extra energy by means of solar radiation during the stirring procedure.
- Do not process any pathogenic materials in open vessels as this may cause the user to inhale vapours. Always work under a suitable extractor hood, even when using closed vessels. Please contact **IKA** customer services if you have any questions regarding work with pathogenic materials or other materials which require closed systems.
- Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Dispersing at high speeds can cause the medium to splash out of the mixing vessel.
- There may be dangerous electrostatic activity between the medium and the dispersing tool.
- The appliance starts up again following a cut in the power supply. Please note that the appliance must be switched off at the on/off switch in this case before the power supply is reconnected.
- The dispersing instrument may not be used in explosive rooms, for mixing dangerous substances or under water.
- Please note that abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on. If there are any questions please contact **IKA**.
- Always disconnect the plug before changing dispersing elements and attaching accessories.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug. The socket for the mains cord must be easily accessible..
- Meter powder not too close to the flange, otherwise powder can be sucked into the drive.
- Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.
- Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.
- The gasket and bearings are made of PTFE and rustproof steel; the following points should therefore be noted: *Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride*

Proper use

The **ULTRA-TURRAX® T 8** is a dispersing unit which is suitable for the production of emulsions and dispersions of minimum quantities in conjunction with the dispersing elements (please order separately) S 8 N - 5 G and S 8 N - 8G. The **ULTRA-TURRAX® T 8** is a hand-held instrument. With the aid of a stand clamp it can also be operated as a stand instrument.

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Unpacking

Please unpack the device carefully and inspect them for any damage. It is important that any damage is detected as soon as the items are unpacked. A record of the relevant facts may be required (railway, mail or carrier).

The device delivery quantity includes: a driving unit, an power supply unit and an operating manual.

Useful information

Because of the high speed of the rotor, the medium to be processed is automatically sucked axially into the dispersing head and then pressed radially through the slots of the rotor-stator-arrangement.


The stator acts at the first approach as the ideal flow breaker. It largely prevents a rotation of the medium and enables an input of big mechanic energies on closest space. The trumpet formation know from stirring is strongly reduced but is not totally avoidable. In the shearing gap between rotor and stator about 1000 times more energy is entered than when stirring the medium.

Setting-up the instrument

Please pay attention to the ambient conditions listed in Technical data (temperature and humidity).

Dispersing elements and drive make an operable dispersing instrument.

- loosen knurled-head screw
- insert dispersing element up to stop
- tighten knurled-head screw

Starting up

Before starting up the instrument, please check the conformity between the voltage indicated on the ID-plate and your mains voltage. If they do not agree, the instrument must not be operated.

- Connect power pack to socket - red LED shines
- Put attaching plug of drive into connector socket of power pack

The attaching plug is equipped with a code pin that ensures the correct polarity. If this code pin is broken, the correct polarity is no longer ensured, and the sense of rotation of the drive may change, which consequently would destroy the dispersing element. In this case the connection lead has to be replaced.

- Put the substance to be processed into e.g. a breaker which has to be secured against torsion by e.g. an IKA-clamp support
- Place the dispersing element into the medium to be processed, and only after that switch on the instrument (at the lowest possible speed)
- Set the requested speed (infinitely variable from 5000 up to 2500 1/min) at the knurling wheel
- In case of solid/liquid mixtures, do not add the solid before the request speed has been reached

Drive unit ULTRA - TURRAX T8

The high-quality drive with a power input of up to max. 100 watts at speeds from 5000 up to 25000 1/min opens new ways in the dispersing technology.

By the infinitely variable speed setting, and by using small dispersing elements with a shaft diameter of 5 and/or 8 mm, smallest quantities can be smoothly and rapidly processed. The drive is equipped with a knurling wheel with a division of 1 to 6 as well as

with a mark on the casing for speed adjustment, so that speeds can be reproducibly set.

Dispersing elements


Dispersing element S 8 N - 5 G

For cleaning purposes, the dispersing element can be dismantled as follows:

- Push rotor shaft (2001) with guide pin upwards out
- draw off bearing bushing (3) from rotor shaft and replace it in case of wear

The assembly is carried out vice-versa.

- | | |
|-------------|------------------------|
| 1 | shaft tube |
| 3 | PTFE - bearing bushing |
| 8 | knurled-head screw |
| 2001 | rotor shaft compl. |


Dispersing tool S 8 N - 8 G

For cleaning purposes, the dispersing element can be dismantled as follows:

- Unscrew rotor (4) from shaft (2)
- push shaft upwards
- draw off bearing bushing (5) from rotor shaft and replace it in case of wear
- push bearing bushing (3) out of shaft and replace it in case of wear

The assembly is carried out vice-versa.

- | | |
|------|------------------------|
| 1 | shaft tube |
| 3 | PTFE - bearing bushing |
| 4 | rotor |
| 5 | bearing bushing |
| 8 | knurled-head screw |
| 2001 | rotor shaft compl. |

Materials contacting the product (both dispersing elements)

- | | |
|-------------------------|---|
| Shaft, rotor and shaft: | anticorrosive steel;
material-No. AISI 316 L |
| Slide-ring bearing: | PTFE |

Maintenance, Cleaning and Sterilisation

Dispersing elements

Cleaning: Immediately after the instrument has been used, the dispersing element must be cleaned so that the threads of rotor and stator do not stick together with remains of the substance. For this purpose the dispersing element is operated in a solvent that removes the remains of the substance, yet does not attack the PTFE-bushings and the anticorrosive steel. Because of the high velocity of flow, rotor and stator are already largely cleaned this way. For a thorough cleaning, however, the dispersing element has to be dismantled, as also the inner walls of the shaft tube are wetted by the medium.

Sterilisation: a) Chemical processing

By bactericidal solvents (formalin, phenol, alcohol) many cases of disinfection can be solved. It is important that the disinfectants are subsequently removed with sterilised water.

b) Germ destruction by moist heat

This is carried out by saturated steam with an excess pressure of 2 bar at 120°C.

c) Germ destruction by hot air

The hot air sterilisation is generally carried out at 160°C and/or 190°C.

The PTFE-bushings (3) and (5) of the dispersing elements have to be supervised. If they are damaged, the medium ascends due to the suction effect of the rotating shaft, and can reach the drive unit this way. So if medium penetrates from the upper lateral bore hole on the shaft tube, any operation has to be interrupted immediately, the PTFE-bushings have to be checked and replaced if required. (Dismantling as described before).

The efficiency of rotor and stator depends on the state of the teeth chamfer. The dispersing effect can be impaired, if due to abrasive media there is a rounding of the teeth.

Drive unit

Maintenance and Cleaning: The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over time.

For cleaning disconnect the main plug.

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.de**, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Accessories

S 8 N - 5 G	Dispersing element
S 8 N - 8 G	Dispersing element
T 8.10	Dispersing station
DTM 12	IKA-TRON Temperature measuring unit
PT 100.23	Temperature measuring sensor
PT 100.24	Temperature measuring sensor

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Technical data


<i>Power pack</i>	nominal voltage	VAC	230 <i>Euro</i>
		VAC	115 <i>USA</i>
	frequency	Hz	50
	output voltage	VDC	12
	current consumption	A	2
	power consumption	VADC	24
	protection class		IP 20
	dimensions WxHxD	mm	84 x 69 x 122
<i>Drive</i>	Weight	kg	1,05
	nominal voltage	VDC	12 ... 18
	power consumption	W	max. 100
	speed	1/min	5000 - 25000
	speed setting		speed setter infinitely variable
	dimensions WxHxD	mm	42 x 44 x 250
	Weight	kg	0,4

<i>S 8 N - 5 G</i>	shaft tube Ø	mm	5
	stator Ø	mm	5
	rotor Ø	mm	3,8
	shearing gapt	mm	0,1
	circumferential speed	m/sec	4,9
	shearing speed	sec ¹	49741
	dept of immersion	mm	85
	use up to	°C	80

<i>S 8 N - 8 G</i>	shaft tube Ø	mm	8
	stator Ø	mm	8
	rotor Ø	mm	6,1
	shearing gapt	mm	0,25
	circumferential speed	m/sec	7,9
	shearing speed	sec ¹	31939
	dept of immersion	mm	110
	use up to	°C	80

List of spare parts drive

Item	Designation
4	clutch
3015	control knob
3016	switch
3018	connecting cable with plug and grommet

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