


INSTRUCTIONS FOR THE *Soap Bubble Flowmeter*

FROM
CHROMATOGRAPHY RESEARCH SUPPLIES

PRODUCT INFORMATION

The CRS Bubble Flowmeter allows accurate calculation of gas flow, by measuring the passage time of a soap film between two marks.



Assembly: Attach rubber tubing to the side-arm of the Bubble Flowmeter. Fill the bulb with soap solution (water and detergent, or Snoop®) and attach to the bottom of the Bubble Flowmeter. Mount the Bubble Flowmeter near the gas supply to be measured, then attach the tubing to the gas supply.

Operation: Squeeze the bulb to raise soap solution above the side arm (gas inlet), then release the bulb to lower soap solution. A film of soap solution will rise in the Bubble Flowmeter. Measure the time of transit of the soap film between two marks, lower and upper. Use a stopwatch for highest accuracy. The two marks give a known volume, and the time of passage allows calculation of volume flow per unit time.

Sample Calculation: Given a 10ml Bubble Flowmeter, with Lower Mark 0.00ml and upper mark 10.0ml; soap film transit time 35 seconds.

$$\frac{10.0\text{ml}}{35\text{sec.}} \times \frac{60\text{sec.}}{1\text{min.}} = 17.1\text{ml/min}$$

REORDER INFORMATION

1ml Soap Bubble Flowmeter	261301	Replacement Bulb, for 25ml	261337
10ml Soap Bubble Flowmeter	261310	Replacement Bulb, for 50ml	261338
25ml Soap Bubble Flowmeter	261325	Replacement Bulb, for 100ml	261339
50ml Soap Bubble Flowmeter	261350	Magnetic Holder (not for 1ml size), pk/2	261333
Triple Stage Bubble Flowmeter	261390	Snoop® Leak Detector	205300
Replacement Bulb, for 1 or 10ml	261335	Stopwatch, Digital	205200




CHROMATOGRAPHY RESEARCH SUPPLIES
2601 TECHNOLOGY DRIVE
LOUISVILLE, KY 40299
PHONE 1-800-327-3800 FAX 502-491-3390
www.chromres.com

INSTRUCTIONS FOR THE *Soap Bubble Flowmeter*

FROM
CHROMATOGRAPHY RESEARCH SUPPLIES

PRODUCT INFORMATION

The CRS Bubble Flowmeter allows accurate calculation of gas flow, by measuring the passage time of a soap film between two marks.



Assembly: Attach rubber tubing to the side-arm of the Bubble Flowmeter. Fill the bulb with soap solution (water and detergent, or Snoop®) and attach to the bottom of the Bubble Flowmeter. Mount the Bubble Flowmeter near the gas supply to be measured, then attach the tubing to the gas supply.

Operation: Squeeze the bulb to raise soap solution above the side arm (gas inlet), then release the bulb to lower soap solution. A film of soap solution will rise in the Bubble Flowmeter. Measure the time of transit of the soap film between two marks, lower and upper. Use a stopwatch for highest accuracy. The two marks give a known volume, and the time of passage allows calculation of volume flow per unit time.

Sample Calculation: Given a 10ml Bubble Flowmeter, with Lower Mark 0.00ml and upper mark 10.0ml; soap film transit time 35 seconds.

$$\frac{10.0\text{ml}}{35\text{sec.}} \times \frac{60\text{sec.}}{1\text{min.}} = 17.1\text{ml/min}$$

REORDER INFORMATION

1ml Soap Bubble Flowmeter	261301	Replacement Bulb, for 25ml	261337
10ml Soap Bubble Flowmeter	261310	Replacement Bulb, for 50ml	261338
25ml Soap Bubble Flowmeter	261325	Replacement Bulb, for 100ml	261339
50ml Soap Bubble Flowmeter	261350	Magnetic Holder (not for 1ml size), pk/2	261333
Triple Stage Bubble Flowmeter	261390	Snoop® Leak Detector	205300
Replacement Bulb, for 1 or 10ml	261335	Stopwatch, Digital	205200



CHROMATOGRAPHY RESEARCH SUPPLIES
2601 TECHNOLOGY DRIVE
LOUISVILLE, KY 40299
PHONE 1-800-327-3800 FAX 502-491-3390
www.chromres.com