# IKA

designed for scientists

**KMO 3 basic**

4

/// ENGLISH

22

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40

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58

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76

/// INDONESIA

94

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ENGLISH

This user guide includes only the main operations and functions.  
You can find and download the **complete operating instructions**  
in your language on our website ([www.ika.com](http://www.ika.com)).

#### Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

## 8

### /// DEVICE SETUP

- › KMO 3 basic . . . . . 8

## 9

### /// SAFETY INSTRUCTIONS

- › About our warning symbols . . . . . 9
- › General information . . . . . 10
- › Device design . . . . . 10
- › Permissible medium / contaminants / side reactions . . . . . 11
- › Procedures during sample runs . . . . . 12
- › Accessories . . . . . 12
- › Power supply / Switching off . . . . . 13
- › About protecting your device . . . . . 14

## 15

### /// PROPER USAGE

- › Scope of delivery . . . . . 15

## 16

### /// START-UP

- › Operation . . . . . 16
- › Operating mode . . . . . 17
- › Automatic reverse rotation . . . . . 18
- › Set-up plate illuminated . . . . . 19
- › Lock . . . . . 19



## 20

### /// TECHNICAL DATA

- › KMO 3 basic in detail . . . . . 20

# Device setup

/// KMO 3 basic

## Explication of symbols

|   |                                   |   |                           |
|---|-----------------------------------|---|---------------------------|
| A | "Power" button                    | H | "+ / -" button            |
| B | Standby LED                       | I | "Reverse rotation" button |
| C | Display                           | J | Threaded hole for stand   |
| D | "Set-up plate illuminated" button | K | Set-up plate              |
| E | "Lock" button                     | L | Power connection          |
| F | "Lock" LED                        | M | Interface USB             |
| G | "Start / Stop" button             |   |                           |

# Safety instructions

/// About our warning symbols

Indicates an imminently hazardous situation, which, if not avoided, **will result in death, serious injury.**

Indicates a potentially hazardous situation, which, if not avoided, **can result in death, serious injury.**

Indicates a potentially hazardous situation, which, if not avoided, **can result in injury.**

Indicates handling, which, if not avoided, **can result in equipment damage.**

**Attention!**  
Indicates the risks of magnetism.

### /// General information

- › **Read the operating manual in its entirety before using the KMO 3 basic and follow the safety instructions.**
- › Keep the operating manual in a location that is easily accessible by everyone in the laboratory.
- › Ensure that only trained staff operates the device.
- › Follow the safety instructions, (laboratory) guidelines, occupational health and safety and accident prevention regulations (OSHA).

#### **Attention – Magnetism!** (📍❤️)

Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).

### /// Device design

#### **Danger!** (⚡)

- › **Do not** operate the device in explosive atmospheres, or under water.
- › Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- › The feet of the device must be clean and undamaged.
- › Check the device and accessories for damage before each use. Do not use damaged components.

#### **Caution!** (🔍)

Device surface is partially made of glass!

- › Glass surface can be damaged by impact.
- › If glass surface is damaged it could cause injury, don't use the device anymore and contact your IKA Service team.

### /// Permissible medium / contaminants / side reactions

#### **Warning!** (!)

Beware of hazards due to:

- › glass breakage.
- › incorrect container size.
- › overfilling of media.
- › unsafe condition of container.
- › Process pathogenic materials only in closed vessels under a suitable fume hood.

Please contact IKA if you have any questions.

#### **Danger!** (⚡)

- › Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- › Please consider any possible contaminations and unwanted chemical reactions.
- › It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- › When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C – 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

*(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)*

### /// Procedures during sample runs

#### **Warning! (!)**

Wear your personal protective equipment in accordance with the hazard category of the media to be processed. There may be a risk from:

- › splashing and evaporation of liquids.
- › ejection of parts.
- › release of toxic or combustible gases.

Reduce speed if:

- › medium splashes out of vessel because the speed is too high.
- › device is not running smoothly.
- › container moves on the working plate.
- › an error occurs.

### /// Accessories

- › Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- › Always disconnect the power before attaching accessories.
- › Observe the operating instructions of the accessories.
- › Accessories must be securely attached to the device so that they cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.

### /// Power supply / Switching off

Place the power supply outside the working area of the device.


- › The voltage stated on the ratings plate must correspond to the mains voltage.
- › The outlet for the power supply must be easily accessible.
- › The device can only be disconnected from main power by unplugging the power supply from the outlet or removing the connector plug at the unit.

#### **Warning! (⚠)**

- › The device will automatically restart in **mode b** following any interruption to the power supply.

### /// About protecting your device

- › The device may only be opened by qualified and IKA approved technicians.
- › Do not cover the device, even partially e.g. with metallic plates or film. This may result in overheating.
- › Protect the device and accessories from bumps and impacts.
- › Keep base plate clean.
- › Observe the minimum distances:  
between devices min. 100 mm,  
between device and wall min. 100 mm,  
above the device min. 800 mm.

## Proper usage

### /// Basics

#### Usage

- › The KMO 3 basic magnetic stirrer has been designed for mixing liquids.

#### Area of use (only indoors)

- › laboratories
- › pharmacies
- › universities
- › schools

#### This device is suitable for use in all areas except

- › residential areas.
- › areas that are connected directly to a low-voltage supply network that also supplies residential areas.

#### The safety of the user cannot be guaranteed

- › if the device is operated with accessories that are not supplied or recommended by the manufacturer.
- › if the device is operated improperly or contrary to the manufacturer's specifications.
- › if the device or the printed circuit board are modified by third parties.

# Start-up

## /// Operation

Check whether the voltage given on the ratings plate corresponds to the available mains voltage. If these conditions have been met, the device is ready for operation when the power supply is plugged in. If these conditions are not met, safe operation is not guaranteed and the device could be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under "Technical data".

### General power supply information

When the power supply adapter is connected to the device, the device changes to the standby mode, the red LED "Standby" **(B)** is on and the display **(C)** will not show any information.

A display self-test will be performed following any interruption to the power supply. The display cycles through the following information:

1. All segments light up.
2. Software version.
3. The selected operating mode.
4. If the drive has not been started the display reads "OFF".
5. If the drive has been started the selected speed is displayed (only in operating mode b).

### Power-up

1. Insert the power supply adapter into the power outlet **(L)**.
2. The device is switched on and off via the "Power" button **(A)**.

### Stirring

1. The stirring function is started and stopped by pressing the "Start / Stop" button **(G)**.
2. Increase or reduce the speed from 0 rpm to 1,500 rpm (in 10 rpm steps) by pressing the "+ / -" buttons **(H)**. Press and hold the "+ / -" button **(H)** may quicken the adjustment of the speed.
  - › *The display **(C)** will only show the set value. If the actual speed is lower than the set speed, the set value on the display will be flashing. The set value on the display will stop flashing once the actual speed reaches the set speed.*

## /// Operating mode

To toggle between operating modes, press the "Reverse rotation" button **(I)** for 2 seconds when the device is powered on.

- › Factory setting: mode A

### Operating mode A

All settings will be stored if the device is disconnected from the power supply. The agitation and "Set-up plate illuminated" functions will be set to "OFF" when the device is powered on. The "Start / Stop" button **(G)** starts the device with the set parameters.

### Operating mode b

All settings will be stored if the device is disconnected from the power supply. The agitation and "Set-up plate illuminated" functions will be set to "ON" or "OFF" when the device is powered on, depending on the previous status of the device.

### /// Automatic reverse rotation

To activate the "Reverse rotation" function, press the "Reverse rotation" button **(I)**. When the "Reverse rotation" function is switched on, a decimal point appears on the LCD display.

The running time (30, 60, 90 or 120 sec) is set by pressing the "Reverse rotation" button **(I)** followed by the "+ / -" buttons **(H)**.

After the preset running time, the drive is switched off and restarted in the opposite direction.

#### Note! (▲)

The activation and deactivation of the „Reverse rotation“ is only possible if the drive is switched off.

### /// Set-up plate illuminated

The illumination brightness can be adjusted by pressing the "Set-up plate illuminated" button **(D)**. After adjustment the brightness value will be shown on the display for 2 seconds.

› Factory setting: Off

| Button pressing time | Set-up plate illuminated | Illumination brightness value * |
|----------------------|--------------------------|---------------------------------|
| 1 x                  | On                       | 50 %                            |
| 2 x                  | On                       | 75 %                            |
| 3 x                  | On                       | 100 %                           |
| 4 x                  | Off                      | 0 %                             |

\* Sequence 50 – 75 – 100 – 0 – 50 – 75 – 100 – 0 – 50 – etc.

### /// Lock

The "Lock" button **(E)** allows the user to lock or unlock the other control buttons (including buttons **D**, **G**, **H** and **I**) of the device.

- › Press the "Lock" button **(E)** for 2 seconds to activate or deactivate the "lock" function.
- › When the "Lock" function is activated, the green LED "Lock" **(F)** is on.
- › When the "Lock" function is deactivated, the green LED "Lock" **(F)** is off.

The device keeps the status of "unlock" or "lock" when the power is switched off and switched on again, which is important for operating mode "b".

# Technical data

/// KMO 3 basic in detail

## Technical data KMO 3 basic

### KMO 3 basic

|  |                    |
|--|--------------------|
| Nominal voltage  | DC 24 / 500 mA     |
| Power input  | 13 W               |
| Power input standby  | 0.8 W              |
| Number of stirring positions   | 1                  |
| Stirring quantity max. (H <sub>2</sub> O)                                | 5 l                |
| Motor rating output  | 9 W                |
| Direction of rotation  | right / left       |
| Speed display set-value  | LED                |
| Speed control  | button             |
| Speed range  | 0 / 10 – 1,500 rpm |
| Setting accuracy speed   | 10 rpm             |
| Speed variation (no load, nominal voltage, at 1,500 rpm, RT: + 25 °C)    | ± 2 %              |
| Stirring bar length  | 20 – 80 mm         |
| Self-heating of the plate by max. stirring (RT: + 22 °C / duration: 1 h) | + 2 K              |
| Set-up plate material  | glass              |
| Set-up plate dimensions  | 140 × 120 mm       |
| Automatic reverse rotation   | yes                |
| Set-up plate illuminated   | yes                |

|  |                     |
|--|---------------------|
| Dimensions (W × H × D)                     | 153 × 62.5 × 227 mm |
| Weight                                     | 1.2 kg              |
| Permissible ambient temperature            | 5 – 40 °C           |
| Permissible relative humidity              | 80 %                |
| Protection class according to DIN EN 60529 | IP 42               |
| USB interface                              | yes                 |

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IKA England LTD.
IKA India Private Limited