

COMMERCIAL CENTER&R-D Center

#203, Building 2, 89 Shengli Rd, Jiangning District, Nanjing 211106, China
Tel.: 0086 25 86216803
Fax: 0086 25 86214463
E-mail: info@citotest.com

PRODUCTION

#339, Beihai Road, Haimen City, Jiangsu 226100, P.R.China
Tel.: 0086 513 82259348
Fax: 0086 513 82110195
E-mail: info@citotest.com.cn

**CITOTEST LABWARE MANUFACTURING CO.,LTD****www.citotest.com**

PRODUCT DATA SHEET

GENERAL INFORMATION

PRODUCT NAME: Microscope Cover Glass

FURTHER DESCRIPTION: Microscope Cover Glass, Circular, Dia.14mm. Thickness 1.0# 0.13-0.16mm. Economical Package.100pcs/box,10x100pcs/inner box,wrapped by aluminum vacumm pack,500box/carton

P/N: 0344-1420

CATEGORY NUMBER: 10210014CE

SPECIFICATION


COLOR: Transparent

MATERIAL: Super white glass

COMPONENT: SiO₂(73%):Na₂O(14%):K₂O(1.5%):CaO(5.82%):MgO(4%):Al₂O₃(1.5%):SO₃(0.15):Fe₂O₃(0.03%)

DIMENSION (mm): Ø14

THICKNESS(mm): 0.13-0.16

PRODUCT PICTURE:**TECHNICAL PROCESS:****INTENDED APPLICATION**

To keep solid specimens pressed flat or liquid samples shaped into a flat layer of even thickness on microscope slides.

PACKING

100pcs/box,10x100pcs/inner box,wrapped by aluminum vacumm pack,500box/carton

BOX:

样品编号 0344-1420 客户编码 4.07.10444 颜色 PANTONE 534C PANTONE 187C 规格尺寸 23X23mm

INNERBOX: [Innerbox Label Size(mm)]



CARTON: [Carton Label Size(mm) 100×130]



STERILIZATION

☐	STERILE	EO	☐	STERILE	R	☐	STERILE	A	☒	
--------------------------	---------	----	--------------------------	---------	---	--------------------------	---------	---	-------------------------------------	--

STORAGE INSTRUCTION

Store in a cool, dry place.

TERM OF VALIDITY

18 months

REMARKS

Conforms with the essential requirements of the European in Vitro Diagnostic 98/79 EC.