



Achat - Mörser mit Pistill

Agate - mortars
with pestle



Maßangaben in mm

Best. Nr.	Øa	Øi	h	t	Vol (ml)
600-M-025	25	20	15	6	2
600-M-035	35	30	18	8	6
600-M-040	40	35	20	10	8
600-M-050	50	40	25	12	10
600-M-060	60	50	25	14	13
600-M-065	65	52	26	15	15
600-M-070	70	58	27	16	20
600-M-075	75	60	30	18	24
600-M-080	80	65	35	20	30
600-M-090	90	75	38	23	40
600-M-100	100	80	40	25	50
600-M-120	120	100	50	30	95
600-M-125	125	105	50	32	100
600-M-150	150	125	65	35	200
600-M-175	175	145	73	45	350
600-M-200	200	165	80	45	550
600-M-110 S	110	90	41	26	65
600-M-115 S	115	93	44	27	80
600-M-130 S	130	113	54	33	125