

Adapter, Joint Reducing & Enlarging
조인트 축소 / 확장 어댑터, LukeGL®

Cat. No.	Model	Inner Socket	Outer Cone	EA Price ₩
A01-83-004	1501214	10/19	14/23	16,600
A01-83-008	1501412	14/23	10/19	18,400
A01-83-016	1501424	14/23	24/40	15,700
A01-83-020	1501429	14/23	29/32	19,800
A01-83-024	1501914	19/26	14/23	17,300
A01-83-028	1501924	19/26	24/40	19,300
A01-83-040	1502414	24/40	14/23	15,700
A01-83-044	1502419	24/40	19/26	19,300
A01-83-048	1502429	24/40	29/32	21,400
A01-83-052	1502434	24/40	34/45	22,500
A01-83-056	1502445	24/40	45/50	34,200
A01-83-060	1502914	29/32	14/23	18,900
A01-83-068	1502924	29/32	24/40	21,400
A01-83-072	1502934	29/32	34/45	27,800
A01-83-076	1502945	29/32	45/50	34,200
A01-83-080	1503424	34/45	24/40	19,300
A01-83-084	1504529	45/50	24/40	34,200



- 요청에 의해 규격 외 사이즈도 공급 가능
- Borosilicate Glass 3.3
- LukeGL®

Glass Reducing Adapter
,Bushing Type
유리 부싱 축소 어댑터

- ASTM E438, Type I, Class A Glass
- Chemglass



Cat. No.	Model	Inner Socket	Outer Cone	EA Price ₩
A01-132-001	1016-27	10/18	14/20	55,500
A01-132-002	1016-28	10/30		55,500
A01-132-901	1016-01	10/30	14/35	77,800
A01-132-905	1016-29	14/20	19/22	54,200
A01-132-909	1016-03	14/35	19/38	64,900
A01-132-004	1016-04	10/30	24/40	77,500
A01-132-005	1016-05	14/35		69,200
A01-132-913	1016-06	19/38	29/42	119,500
A01-132-917	1016-08	14/35		118,900
A01-132-921	1016-09	19/38	115,200	
A01-132-925	1016-10	24/40	34/45	74,500
A01-132-029	1016-12			96,700
A01-132-933	1016-13	29/42	119,200	
A01-132-937	1016-14	24/40	40/50	159,300
A01-132-941	1016-17		164,200	
A01-132-945	1016-18	29/42	45/50	156,000
A01-132-949	1016-19	34/45		154,900

PTFE Reducing Adapter
,Bushing Type
테프론 부싱 축소 어댑터

- 재질 : PTFE, -200 ~ 280°C 내열
- PTFE 테프론 재질로 최상의 내열, 내화학성을 가짐
- Cowie



Full Length

Cat. No.	Model	Inner Socket	Outer Cone	EA Price ₩
A01-450-100	021.01014	10/30	14/35	16,600
A01-450-102	021.01019		19/38	19,100
A01-450-104	021.01419	14/35	19/38	21,400
A01-450-106	021.01424		24/40	28,100
A01-450-109	021.01924	19/38	24/40	28,100
A01-450-111	021.01929		29/42	33,600
A01-450-113	021.01934	24/40	34/45	41,400
A01-450-115	021.02429		29/42	39,800
A01-450-117	021.02434	29/42	34/45	43,700
A01-450-119	021.02445		45/50	58,000
A01-450-121	021.02934	29/42	34/45	41,300
A01-450-123	021.02945		45/50	54,900
A01-450-125	021.03445	34/45	45/50	65,500

Medium Length

Cat. No.	Model	Inner Socket	Outer Cone	EA Price ₩
A01-450-131	021.1014	10/19	14/29	17,500
A01-450-134	021.1019		19/26	19,400
A01-450-137	021.1419	14/23	19/26	20,800
A01-450-140	021.1424		24/29	25,900
A01-450-143	021.1924	19/26	24/29	25,900
A01-450-147	021.1929		29/32	32,600
A01-450-150	021.1934	24/29	34/45	36,700
A01-450-153	021.2429		29/32	32,600
A01-450-156	021.2434	29/32	34/45	38,900
A01-450-159	021.2445		45/40	43,500
A01-450-162	021.2934	29/32	34/45	36,700
A01-450-165	021.2945		45/40	41,200
A01-450-168	021.3445	34/45	45/40	55,400

PTFE Ground Joint Reducing Set
테프론 조인트 어댑터 세트

Cat. No.	Model	Description	EA Price ₩
A01-460-855	H981-14	NS14, NS19, NS24, NS29, NS34, NS45, NS60	195,100



- 재질 : PTFE, -200 °C to 250°C 내화학 및 내열
- For bridging different cone and socket size e.g. a NS 14 cone can functionally be placed into a NS 45 socket
- BOLA