

I Absorption Cell, 2-Side Polished, Economy Type

경제형 흡광 셀, 2면 투명

• YJ Optical

Standard Cell with PTFE Cover, Type 1

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩	Optical Glass (340~2500 nm)		EA Price ₩
		Cat No.	Model		Cat No.	Model	
1 mm	0.4 ml	S18-203-002	1UV-1	69,400	S18-203-802	1GL-1	14,900
2 mm	0.7 ml	S18-203-004	1UV-2	56,400	S18-203-804	1GL-2	12,700
5 mm	1.7 ml	S18-203-006	1UV-5	45,500	S18-203-806	1GL-5	12,800
10 mm	3.5 ml	S18-203-008	1UV-10	45,500	S18-203-808	1GL-10	7,800
20 mm	7.0 ml	S18-203-010	1UV-20	56,400	S18-203-810	1GL-20	13,000
40 mm	14.0 ml	S18-203-014	1UV-40	84,600	S18-203-814	1GL-40	15,100
50 mm	17.5 ml	S18-203-016	1UV-50	97,600	S18-203-816	1GL-50	21,000
100 mm	35.0 ml	S18-203-018	1UV-100	225,500	S18-203-818	1GL-100	40,000

• Inside Width : 10 mm

• Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm

Standard Cell with PTFE Stopper, Type 21

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩	Optical Glass (340~2500 nm)		EA Price ₩
		Cat No.	Model		Cat No.	Model	
1 mm	0.4 ml	S18-203-020	21UV-1	147,500	S18-203-820	21GL-1	116,000
2 mm	0.7 ml	S18-203-022	21UV-2	147,500	S18-203-822	21GL-2	116,000
5 mm	1.7 ml	S18-203-024	21UV-5	65,100	S18-203-824	21GL-5	30,000
10 mm	3.5 ml	S18-203-026	21UV-10	65,100	S18-203-826	21GL-10	26,400
20 mm	7.0 ml	S18-203-027	21UV-20	75,900	S18-203-827	21GL-20	29,800
40 mm	14.0 ml	S18-203-028	21UV-40	101,900	S18-203-828	21GL-40	46,000
50 mm	17.0 ml	S18-203-030	21UV-50	117,100	S18-203-830	21GL-50	50,000

• Inside Width : 10 mm

• Outside Size : w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm

Semi-Micro Cell

Cat. No.	Model	Type	Description	EA Price ₩
S18-203-050	9UV-10	9	Semi-Micro Cell with PTFE Cover	65,100
S18-203-062	9MUV-10	9M	Black Semi-Micro Cell with PTFE Cover	60,700
S18-203-070	29UV-10	29	Semi-Micro Cell with PTFE Stopper	73,700
S18-203-082	29MUV-10	29M	Black Semi-Micro Cell with PTFE Stopper Cover	67,200

• UV Quartz (190 ~ 2500 nm)

• Inside Width : 4 mm (9UV-10, 29UV-10) / 10 mm (9MUV-10, 29MUV-10)

• Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm (9UV-10, 9MUV-10)
w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm (29UV-10, 29MUV-10)

• Volume : 1.4 ml, Light Path : 10 mm

Micro Cell

Cat. No.	Model	Type	Description	EA Price ₩
S18-203-116	18UV-10	18	Micro Cell with PTFE Cover	65,100
S18-203-128	18MUV-10	18M	Black Micro Cell with PTFE Cover	60,700
S18-203-140	30UV-10	30	Micro Cell with PTFE Cove	73,700
S18-203-152	30MUV-10	30M	Black Micro Cell with PTFE Cover	67,200

• UV Quartz (190 ~ 2500 nm)

• Inside Width : 2 mm

• Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm (18UV-10, 18MUV-10)
w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm (30UV-10, 30MUV-10)

• Volume : 0.7 ml, Light Path : 10 mm



Type 1



Type 21



Type 9



Type 9M



Type 29



Type 29M



Type 18



Type 18M



Type 30



Type 30M