


Absorption Cell, 2-Side Polished, Economy Type 에서 계속...

Rectangular Absorption meter Cell, Type 101

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩
		Cat No.	Model	
2.5 mm	2.2 ml	S18-203-214	101UV-2.5	73,700
5 mm	4.4 ml	S18-203-216	101UV-5	58,600
10 mm	8.8 ml	S18-203-218	101UV-10	58,600
20 mm	17.5 ml	S18-203-220	101UV-20	75,900
40 mm	35.0 ml	S18-203-222	101UV-40	106,300
50 mm	43.8 ml	S18-203-224	101UV-50	117,100
100 mm	87.5 ml	S18-203-226	101UV-100	264,600
200 mm	175.0 ml	S18-203-228	101UV-200	529,200

- Inside Width : 2 mm
- Outside Size : w28×h40 mm, d (깊이) = Light Path + 4 mm



Type 101

Cylindrical Cell with PTFE Stopper, Type 32

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩
		Cat No.	Model	
10 mm	2.8 ml	S18-203-236	32UV-10	52,000
20 mm	5.6 ml	S18-203-238	32UV-20	67,200
50 mm	14.0 ml	S18-203-240	32UV-50	93,300
100 mm	28.0 ml	S18-203-242	32UV-100	132,300

- Inside Size : Ø19 mm
- Outside Size : Ø22×h30 mm, d (깊이) = Light Path + 2.5 mm
- 모두 1개의 Stopper 만을 가짐



Type 32

Fluorescence Cell, 4-Side Polished

형광 셀, 4면 투명

- YJ Optical

Standard Fluorometer Cell with PTFE Cover, Type 3

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩	Optical Glass (340~2500 nm)		EA Price ₩
		Cat No.	Model		Cat No.	Model	
10 mm	3.5 ml	S18-203-336	3UV-10	124,000	S18-203-972	3GL-10	26,400
20 mm	7.0 ml	S18-203-338	3UV-20	99,800	-	-	-

- Inside Width : 10 mm
- Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm



Type 3



Type 23

Standard Fluorometer Cell with PTFE Stopper, Type 23

Light Path	Volume	UV Quartz (190~2500 nm)		EA Price ₩	Optical Glass (340~2500 nm)		EA Price ₩
		Cat No.	Model		Cat No.	Model	
10 mm	3.5 ml	S18-203-348	23UV-10	93,300	S18-203-984	23GL-10	36,400

- Inside Width : 10 mm
- Outside Size : w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm



Type 517



Type 29F

Semi-Micro Fluorometer Cell

Cat. No.	Model	Description	EA Price ₩
S18-203-362	517UV-10	with PTFE Cover, Type 517	117,100
S18-203-364	29FUV-10	with PTFE Stopper, Type 29F	132,300

- UV Quartz (190 ~ 2500 nm)
- Inside Width : 4 mm
- Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm (517UV-10)
w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm (29FUV-10)
- Volume : 1.4 ml, Light Path : 10 mm

Micro Fluorometer Cell

Cat. No.	Model	Description	EA Price ₩
S18-203-366	18FUV-10	with PTFE Cover, Type 18F	117,100
S18-203-368	30FUV-10	with PTFE Stopper, Type 30F	132,300

- UV Quartz (190~2500 nm)
- Inside Width : 2 mm
- Outside Size : w12.5×h45 mm, d (깊이) = Light Path + 2.5 mm (18FUV-10)
w12.5×h49 mm, d (깊이) = Light Path + 2.5 mm (30FUV-10)
- Volume : 1.4 ml, Light Path : 10 mm



Type 18F



Type 30F