

Absorption Ultra-Micro Cell 울트라 마이크로 흡광 셀

• Hellma

Type : 105.200-QS, Quartz

Cat. No.	Model	Path	Height	Chamber	용량	EA Price ₩
S18-69-498	105.200-85-QS	10mm	8.5mm	160 ul	180 ul	861,600
S18-69-499	105.200-15-QS		15.0mm			861,600

• Aperture : w2×h8 mm

Type : 105.201-QS, Quartz

Cat. No.	Model	Path	Height	Chamber	용량	EA Price ₩
S18-69-455	105.201-85-QS	10mm	8.5mm	100 ul	120 ul	861,600
S18-69-459	105.201-15-QS		15.0mm			861,600

• Aperture : w2×h5 mm

Type : 105.202-QS, Quartz

Cat. No.	Model	Path	Height	Chamber	용량	EA Price ₩
S18-69-466	105.202-85-QS	10mm	8.5mm	50 ul	70 ul	861,600
S18-69-469	105.202-15-QS		15.0mm			861,600

• Aperture : w2×h2.5 mm

Type : 105.210-QS, Quartz

Cat. No.	Model	Path	Height	Chamber	용량	EA Price ₩
S18-69-492	105.210-1085-QS	10mm	8.5mm	5 ul	10 ul	1,459,800
S18-69-496	105.210-1015-QS		15.0mm			1,459,800

• Aperture : Ø0.8 mm



105-200-15-QS



105-201-85-QS



105-202-15-QS



105-210-1015-QS

Absorption Cell for Magnetic Stirrer 자석교반용 흡광 셀

• Hellma

Type : 109.000-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-502	109.000-10-QS	10 mm	3.5 ml	518,800
S18-69-506	109.004-10-QS		1.5 ml	581,600

• Base Thickness : 5 mm

• With PTFE Lid

Type : 119.000-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-510	119.000-10-QS	10 mm	3.5 ml	631,800
S18-69-514	119.004-10-QS	10 ml	1.5 ml	694,600

• Base Thickness : 5 mm

• With PTFE Stopper



109.000-10-QS



119.000-10-QS

Absorption Sealable Cell 스크류 캡 흡광 셀

• Hellma

Type : 117.100-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-518	117.100-10-QS	10 mm	3.5 ml	692,000

• Base Thickness : 1.5 mm, GL 14

Type : 117.104-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-522	117.104-10-QS	10 mm	1.4 ml	807,600

• Base Thickness : 1.25 mm, GL 14



117.100-10-QS



117.104-10-QS