

## Absorption Flow-Through Cell 연속식 흡광 셀

• Hellma

Type : 130-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	Aperture	용량	EA Price ₩
S18-69-560	130-10-QS	10 mm	9.5×33 mm	3.2 ml	802,500

• Size : h45×w12.5×d12.5 mm

Type : 137-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	Aperture	용량	EA Price ₩
S18-69-564	137-10-QS	10 mm	9×20 mm	2.6 ml	930,800

• Size : h45×w12.5×d12.5 mm

Type : 170-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	Aperture	용량	EA Price ₩
S18-69-576	170-1-QS	1 mm	3.5×17.5 mm	120 ul	1,273,900
S18-69-580	170-2-QS	2 mm		240 ul	1,273,900

• Size : h35×w12.5×d12.5 mm

Type : 175.000-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	Center Height	용량	EA Price ₩
S18-69-588	175-15-10-QS	10 mm	15 mm	750 ul	975,900

• Size : h38.5×w12.5×d12.5 mm (175-85-10-QS)

h45×w12.5×d12.5 mm (175-15-10-QS)

Type : 170.700-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-592	170.700-0.1-QS	0.1 mm	6 ul	1,857,300

• Size : w12.5×d12.5×h35 mm



130-10-QS



137-10-QS



170-1-QS

## Absorption Cylindrical Cell 실린더형 흡광 셀, Type 120-QS

Cat. No.	Model	Path	용량	EA Price ₩
S18-69-530	120-1-QS	1 mm	280 ul	609,200
S18-69-534	120-2-QS	2 mm	560 ul	604,200
S18-69-542	120-10-QS	10 mm	2.8 ml	355,900
S18-69-546	120-20-QS	20 mm	5.6 ml	398,300
S18-69-550	120-50-QS	50 mm	14.0 ml	579,000

• With PTFE Stopper (#120-50-QS는 2 PTFE Stopper임)

• Operating Range : 200 ~ 2500 nm

• Size : OD22 mm, ID19 mm • Hellma



120-50-QS

## Secondary Spectrophotometer Calibration Standard 2차 흡광 셀 측정 표준

Cat. No.	Model	Description	EA Price ₩
S18-69-755	666-000	Calibration Standards	3,479,900

• Wavelength와 Photometric Accuracy의 확인용

• Absorbance in visible range: 440, 465, 546, 590, 635nm wavelength

• Wavelength in the UV and Visible Range: 279.3, 360.9, 453.5, 536.3, 637.6nm

• Hellma

### 구성 Component

Cat. No.	Model	재질	Nominal Absorban	EA Price ₩
S18-69-911	666-F1	Holmium Oxide Glass Filter	-	937,800
S18-69-915	666-F2	Grey Glass Filter NG 11	0.25	937,800
S18-69-919	666-F3	Grey Glass Filter NG 5	0.5	937,800
S18-69-923	666-F4	Grey Glass Filter NG 4	1.0	937,800
S18-69-927	666-F0	Empty Mount	-	222,800

