

Fluorescence Macro Cell

표준 형광 셀

- Base Thickness : 1.25 mm
- Hellma

Type : 101-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-644	101-10-QS	10×10 mm	3.5 ml	410,800

- With PTFE Cover

Type : 111-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-654	111-10-QS	10×10 mm	3.5 ml	546,500

- With PTFE Stopper



101-10-QS



111-10-QS

Fluorescence Semi-Micro Cell

세미 마이크로 형광 셀

- Base Thickness : 1.25 mm
- Hellma

Type : 104F-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-664	104F-10-QS	10×4 mm	1.4 ml	749,900

- With PTFE Cover

Type : 114F-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-674	114F-10-QS	10×4 mm	1.4 ml	875,900
S18-100-011	7613402	PTFE Stopper for cuvettes		11,700

- With PTFE Stopper



104F-10-QS



114F-10-QS

Fluorecne Micro Cell

마이크로 형광셀

- Hellma

Type : 104.002F-QS, Quartz
Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-684	104.002F-10-QS	10×2 mm	700 ul	835,200

- Base Thickness : 1.25 mm
- With PTFE Cover

Type : 108.002F-QS, Quartz
Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-694	108.002F-10-QS	10×2 mm	500 ul	870,400

- Base Thickness : 9 mm
- With PTFE Cover

Type : 115F-QS, Quartz / Operating Range : 200 ~ 2500 nm

Cat. No.	Model	Light Path	용량	EA Price ₩
S18-69-698	115F-10-QS	10×2 mm	400 ul	898,000

- Base Thickness : 1.25 mm
- With PTFE Stopper



104.002F-10-QS



108.002F-10-QS



115F-10-QS