

Disposable UV Micro Cuvette UV-VIS용 일회용 마이크로 큐벳

- Specially Designed UV-Vis Micro Cuvette : 70 ~ 850 ul
- For Photometric Determination of Protein, ssDNA, dsDNA, RNA, Oligonucleotide
- Wavelength Range : 220 ~ 900 nm
- Size : w/d12.5×h45 mm
- 2면 투명
- Brand GmbH

Non-Sterile

Cat. No.	Model	용량	Center Height	Window	PK Price ₩	
S19-59-111	759200	70 ~ 850 ul	8.5 mm ①	2면 투명	100ea	62,300
S19-59-115	759220		15 mm ②		100ea	58,200

Sterile, Single Wrapped


Cat. No.	Model	용량	Center Height	Window	PK Price ₩	
S19-59-121	759215	70 ~ 850 ul	8.5 mm	2면 투명	100ea	120,700
S19-59-125	759235		15 mm		100ea	120,700

- Free of DNase, DNA and RNase

Cap for UV Micro Cuvette ③

Cat. No.	Model	Type	재질	색상	PK Price ₩	
S19-59-131	759240	Round	PE	Blue	100ea	27,700
S19-59-135	759241			Yellow	100ea	27,700

- UV Micro Cuvette과 같이 -20°C 까지 Sample 보관 가능



759200

759220



759240

Disposable Standard Cuvette 일회용 큐벳

- 2 Type : PS 재질, PMMA 재질
 - PS Wavelength Range : 340 ~ 900 nm
 - PMMA Wavelength Range : 300 ~ 900 nm
- Light Path : 10 mm
- Size : w/d12.5×h45 mm
- Brand GmbH

PS 재질, Operating Range : 340 ~ 900 nm

Cat. No.	Model	용량	Type	Window	PK Price ₩	
S18-59-010	759005	2.5 ~ 4.5 ml	Macro	2면 투명	100ea	19,900
S18-59-014	759015	1.5 ~ 3.0 ml	Semi-Micro		100ea	21,000
S18-59-016	759030	2.5 ~ 4.5 ml	Macro	4면 투명	100ea	30,800


PMMA 재질, Operating Range : 300 ~ 900 nm

Cat. No.	Model	용량	Type	Window	PK Price ₩	
S18-59-018	759105	2.5 ~ 4.5 ml	Macro	2면 투명	100ea	28,700
S18-59-022	759115	1.5 ~ 3.0 ml	Semi-Micro		100ea	30,000

Disposable UV Cuvette UV용 일회용 큐벳, 230~900 nm

Cat. No.	Model	용량	Type	Window	PK Price ₩	
S18-59-026	759150	1.5 ~ 3.0 ml	Semi-Micro	2면 투명	100ea	80,000
S18-59-030	759170	2.5 ~ 4.5 ml	Macro		100ea	72,500
S18-59-034	759125	2.5 ~ 4.5 ml	Macro	4면 투명	100ea	83,500

- Water Analysis, Chemistry, Life Science Application에 사용하기에 알맞음
- 대부분의 솔벤트, 산, 알칼리 용액 사용 가능
- 석영 재질의 Cuvette에 비해 내용물 오염이 적음
- Light Path : 10 mm
- Wavelength Range : 230 ~ 900 nm
- Size : w/d12.5×h45 mm
- Brand GmbH



759005, 759105, 759170

759015, 759115, 759150



759030, 759125