

Disposable Tricorn Beaker 일회용 PP 비이커

Cat. No.	Model	용량	Size (상/하 $\varnothing \times h$)	CS Price ₩	
B05-53-008	DTB50	50 ml	47/32×57 mm	100ea	38,600
B05-53-016	DTB100	100 ml	59/39×71 mm	100ea	51,800
B05-53-024	DTB250	250 ml	74/50×88 mm	100ea	73,500
B05-53-032	DTB400	400 ml	84/52×114 mm	100ea	80,700
B05-53-040	DTB800	800 ml	106/65×137 mm	100ea	112,000
B05-53-048	DTB1000	1 L	115/72×143 mm	100ea	154,200

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- 반투명 내화학 재질
- 일회성의 보급형 비이커이나 재사용 가능
- 3개의 Pouring Spout와 몰드 눈금
- LabPlasti®



Disposable Flexible Beaker 일회용 PP 비이커

Cat. No.	Model	용량	Size	CS Price ₩	
B05-24-916	DFB100	100 ml	$\varnothing 54/47 \times h72$ mm	1000ea	267,700
B05-24-924	DFB250	200 ml	$\varnothing 73/63 \times h84$ mm	1000ea	535,400
B05-24-932	DFB400	300 ml	$\varnothing 84/67 \times h99$ mm	500ea	348,000
B05-24-940	DFB800	500 ml	$\varnothing 99/83 \times h121$ mm	250ea	200,700
B05-24-948	DFB1000	1 L	$\varnothing 122/104 \times h145$ mm	100ea	139,200
B05-24-952	DFB2000	2 L	$\varnothing 142/124 \times h197$ mm	100ea	321,200

- PP 재질의 Flexible 비이커
- Universal



Disposable Polystyrene Beaker 일회용 플라스틱 비이커

Cat. No.	Model	용량	Size (상/하 $\varnothing \times h$)	PK Price ₩	CS Price ₩
B05-134-010	235336	5 ml	33/19×19 mm	100ea 24,900	1000ea 236,500
B05-134-014	235337	10 ml	36/22×25 mm	100ea 23,100	1000ea 215,700
B05-134-018	235338	20 ml	44/28×31 mm	100ea 51,800	500ea 246,000
B05-134-022	235339	50 ml	60/35×50 mm	100ea 80,200	500ea 380,900
B05-134-024	235340	150 ml	84/54×63 mm	50ea 48,700	100ea 91,000
B05-134-028	235341	250 ml	92/65×76 mm	50ea 60,500	100ea 113,100

- 재질 : PS, -20 ~ 70°C 내열, 반투명
- 4방향의 Spout을 가지고 있음
- 5, 10, 20, 50 ml는 Atomic absorption, Flame photometers, Blood cell counter 용도에 알맞음
- Qorpak



ETFE Beaker ETFE 테프론 비이커

Cat. No.	Model	용량	Size	EA Price ₩
B05-59-712	90224	25 ml	$\varnothing 32 \times h50$ mm	16,900
B05-59-716	90228	50 ml	$\varnothing 39 \times h59$ mm	22,300
B05-59-720	90238	100 ml	$\varnothing 50 \times h72$ mm	30,800
B05-59-724	90248	250 ml	$\varnothing 67 \times h96$ mm	52,300
B05-59-728	90252	400 ml	$\varnothing 77 \times h109$ mm	72,500
B05-59-732	90254	500 ml	$\varnothing 88 \times h122$ mm	84,800
B05-59-736	90256	600 ml	$\varnothing 91 \times h125$ mm	98,200
B05-59-740	90262	1000 ml	$\varnothing 105 \times h143$ mm	139,900

- 재질 : ETFE, -100 ~ 150°C 내열
- ETFE 테프론 재질로 뛰어난 내화학성을 가짐
- Brand GmbH

