


Cell Compartment Spacer 셀 스페이서

Cat. No.	Model	Description	EA Price ₩
S18-77-611	SPA11	Spacer for 1 mm Light Path Cell	84,700
S18-77-615	SPA12	Spacer for 2 mm Light Path Cell	84,700
S18-77-619	SPA15	Spacer for 5 mm Light Path Cell	84,700

- 재질 : Aluminum Black Anodising Spacer
- 용도 : Standard 10 mm Cuvette Holder에 Light Path 1,2,5 mm Cell을 장착하기 위한 Spacer(Type 1과 Type 2)
- Starna



Cell Mount 셀 마운트

Cat. No.	Model	적용 셀	EA Price ₩
S18-203-A20	CH/2049	Light path 3mm 이하 사각 Cell용	131,800
S18-203-A23	A23	사이즈 7×7 mm 사각 Cell 용	328,900

- Starna : CH/2049
- ICL : A23



CH/2049



A23

Square Cuvette Holder 사각 큐벳 홀더

Cat. No.	Model	적용 셀	Aperture	EA Price ₩
S18-87-209	0503	10 mm Light path용 Cell	20×10 mm	388,400
S18-87-216	2466	1-50 mm Light path용 Cell, Adjustable		717,300

- 모든 IR 또는 FTIR Spectrophotometer들의 2"×3" Universal Slid Amount
- ICL



0503




2466

Cylinder Cell Holder 실린더형 큐벳 홀더

Cat. No.	Model	적용 큐벳	Aperture	EA Price ₩
S18-87-305	0023-4607	10 mm Light path용 Cell	Ø19 mm	345,800
S18-87-312	0023-2984	50 mm Light path용 Cell		337,300
S18-87-319	0023-502	100 mm Light path용 Cell		324,600

- ICL



Cuvetter Cleaner, Cell Washer 큐벳 세척기

Cat. No.	Model	Position	Description	EA Price ₩
W07-10-552	P65S	1	Single Cell Washer	536,100
W07-10-560	P65D	2	Double Cell Washer	315,200
F04-985-101	412022500		Filtering Flask, 500 ml	61,300

- 큐벳의 세척이나 건조용
- OD 20 mm 이내의 작은 유리류의 세척용
- 구성 : Stopper 포함
- 500 ml Filtering Flask는 별매
(Size : Ø105×h185 mm, Neck ID : 35 mm)
- NSG

