


Colored Micro Magnetic Stir Bar 칼라 마이크로 마그네틱 바

- PTFE Covered
- 내열온도 : -200 ~ 280°C
- Cowie



Ø3 x L3 mm Micro Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-701	001.803B	Blue	2,200
M01-450-703	001.803R	Red	2,200
M01-450-705	001.803Y	Yellow	2,200

Ø3 x L6 mm Micro Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-707	001.806B	Blue	1,300
M01-450-711	001.806R	Red	1,300
M01-450-715	001.806Y	Yellow	1,300

Ø3 x L10 mm Micro Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-727	001.810B	Blue	1,100
M01-450-731	001.810R	Red	1,100
M01-450-735	001.810Y	Yellow	1,100

Ø8 x L13 mm Micro Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-901	001.513B	Blue	2,200
M01-450-903	001.513R	Red	2,200
M01-450-905	001.513Y	Yellow	2,200

Assorted Magnetic Stir Bar Set 마그네틱 바 셋

- 플라스틱 Case 포장
- Cowie



Cat. No.	Model	규격	EA Price ₩
M01-450-060	CW-P20	Plan type Set, 20ea	83,300
M01-450-070	CW-O18	Octagonal Set, 16ea	57,200
M01-450-080	CW-E14	Egg Set, 14ea	61,000
M01-450-090	CW-A52	Assorted Set, 50ea	201,500

- CW-P20 : 4.5×15, 7×20, 7×30, 7×35, 8×25, 8×40, 8×50, 10×60, 10×70, 10×80 mm (각 2개)
- CW-O18 : 8×13, 8×15, 8×22, 8×25, 8×28, 8×41, 8×51, 10×64 mm (각 2개)
- CW-E14 : 5×10, 6×15, 10×20, 10×25, 16×30, 16×35, 20×40 mm (각 2개)
- CW-A52 : CW-P20 + CW-O18 + CW-E14

Colored Magnetic Stir Bar Octagonal Shaped 칼라 8변형 마그네틱 바

- PTFE Covered
- 내열온도 : -200 ~ 280°C
- Cowie



Ø8 x L15 mm Octagonal Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-915	001.515B	Blue	3,100
M01-450-916	001.515R	Red	3,100
M01-450-917	001.515Y	Yellow	3,100

Ø8 x L25 mm Octagonal Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-709	001.525B	Blue	3,300
M01-450-713	001.525R	Red	3,300
M01-450-717	001.525Y	Yellow	3,300

Ø8 x L38 mm Octagonal Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-729	001.538B	Blue	4,600
M01-450-733	001.538R	Red	4,600
M01-450-737	001.538Y	Yellow	4,600

Ø8 x L51 mm Octagonal Bar


Cat. No.	Model	Color	EA Price ₩
M01-450-749	001.551B	Blue	5,300
M01-450-753	001.551R	Red	5,300
M01-450-757	001.551Y	Yellow	5,300

Ø13 x L75 mm Octagonal Bar

Cat. No.	Model	Color	EA Price ₩
M01-450-769	001.575B	Blue	11,100
M01-450-773	001.575R	Red	11,100
M01-450-777	001.575Y	Yellow	11,100

Magnetic Retriever 마그네틱 리트리바

- Cowie



PTFE Retriever

Cat. No.	Model	길이	EA Price ₩
M01-450-001	004.150	150 mm	19,900
M01-450-005	004.250	250 mm	24,400
M01-450-009	004.350	350 mm	27,700
M01-450-013	004.450	450 mm	32,000

PP Retriever

Cat. No.	Model	길이	EA Price ₩
M01-450-033	004.1300	300 mm	12,000
M01-450-037	004.1350	350 mm	12,400
M01-450-041	004.1450	450 mm	13,100