

Turbo Cylindrical Shaped Magnetic Stir Bar 원통형 터보 마그네틱 바

- PTFE Coated (코발트 자석)
- 내열온도 : -200 ~ 280°C
- Cowie



Turbo
Cylindrical Shaped

Cat. No.	Model	Size	EA Price ₩
M07-450-001	001.108.RE	∅3.0×L8 mm	3,500
M07-450-005	001.112.RE	∅4.5×L12 mm	4,300
M07-450-009	001.120.RE	∅6.0×L20 mm	5,400
M07-450-013	001.125.RE	∅6.0×L25 mm	7,400
M07-450-017	001.130.RE	∅6.0×L30 mm	7,900
M07-450-021	001.140.RE	∅8.0×L40 mm	9,300
M07-450-025	001.150.RE	∅8.0×L50 mm	9,700
M07-450-029	001.160.RE	∅10.0×L60 mm	12,200

Turbo Octagonal Shaped Magnetic Stir Bar 팔면체형 터보 마그네틱 바

- PTFE Coated (코발트 자석)
- 내열온도 : -200 ~ 280°C
- Cowie




Turbo
Octagonal Shaped

Cat. No.	Model	Size	EA Price ₩
M07-450-222	001.513.3.RE	∅3.0×L13 mm	3,700
M07-450-226	001.513.RE	∅8.0×L13 mm	5,400
M07-450-230	001.515.RE	∅8.0×L15 mm	5,800
M07-450-234	001.525.RE	∅8.0×L25 mm	7,700
M07-450-238	001.538.RE	∅8.0×L38 mm	7,900
M07-450-242	001.551.RE	∅8.0×L51 mm	9,100
M07-450-246	001.564.8.RE	∅8.0×L64 mm	10,500

Turbo Magnetic Stir Bar Set 터보 마그네틱 바 셋

- 플라스틱 Case 포장
- Cowie




Stir Bar Set

Cat. No.	Model	Size	EA Price ₩
M01-450-100	CW-TC14	Turbo Cylindrical Set, 14ea	97,800
M01-450-110	CW-TE14	Turbo Elliptical Set, 10ea	82,000

- CW-TC14 : ∅8×L50, ∅10×L60 (각 1개), ∅3×L8, ∅4.5×L12, ∅6×L20, ∅6×L25, ∅6×L30, ∅8×L40 mm (각 2개)
- CW-TE14 : ∅24×L50 (각 1개), ∅6×L10, ∅10×L15, ∅14×L25 mm (각 3개)

Turbo Elliptical Shaped Magnetic Stir Bar 타원형 터보 마그네틱 바

- PTFE Coated (코발트 자석)
- 내열온도 : -200 ~ 280°C
- Cowie



Turbo
Elliptical Shaped

Cat. No.	Model	Size	EA Price ₩
M07-450-150	001.2610.RE	∅6.0×L10 mm	4,100
M07-450-154	001.2615.RE	∅10.0×L15 mm	6,400
M07-450-158	001.2625.RE	∅14.0×L25 mm	8,600
M07-450-162	001.2650.RE	∅24.0×L50 mm	23,700
M07-450-166	001.2670.RE	∅28.0×L70 mm	43,300

Turbo Egg Shaped Magnetic Stir Bar 팔면체형 마그네틱 바

- PTFE Coated (코발트 자석)
- 내열온도 : -200 ~ 280°C
- Cowie



Turbo Egg Shaped

Cat. No.	Model	Size	EA Price ₩
M07-450-301	001.610.RE	∅5.0×L10 mm	3,200
M07-450-305	001.615.RE	∅6.0×L15 mm	4,100
M07-450-309	001.620.RE	∅10.0×L20 mm	8,000
M07-450-313	001.625.RE	∅12.0×L25 mm	7,200
M07-450-317	001.635.RE	∅16.0×L35 mm	11,400
M07-450-321	001.640.RE	∅20.0×L40 mm	16,300
M07-450-325	001.650.RE	∅20.0×L50 mm	26,400
M07-450-329	001.664.RE	∅20.0×L64 mm	20,700

Turbo Magnetic Retriever 마그네틱 터보 리트리버

- PTFE Coated (코발트 자석)
- 내열온도 : -200 ~ 280°C
- Cowie



Turbo Retriever

Cat. No.	Model	Size	EA Price ₩
M07-450-560	004.150.RE	∅10.0×L150 mm	21,200
M07-450-564	004.250.RE	∅10.0×L250 mm	26,600
M07-450-568	004.350.RE	∅10.0×L350 mm	29,300
M07-450-572	004.450.RE	∅10.0×L450 mm	32,100