

Magnetic Stirrer, MIX 1 eco 자력 교반기

Cat.No.	Model	Speed Range	Capacity	EA Price ₩
S08-75-070	30101	120 ~ 1400 rpm	1 ~ 3000 ml	812,800

- Suitable for the operation with aggressive media
- Ultra Flat, Wear-free 자력 교반기, All Stainless Steel 재질
- Start가 부드러우며 가속시 바의 틈 현상이 없음 • 2 mag

Cat. No.	S08-75-070
Model	30101
Stirring Point	1
Stirring Volume	1 ~ 3000 ml
Speed Range	120 ~ 1400 rpm
Stirring Power	10 W
Material Housing	Stainless Steel
Operation Conditions	-10 up to + 50°C (at 80% Humidity)
Electrical Data	100 ~ 240 V / 50 ~ 60 Hz / 0.7 A
Size / Weight	w120×d120×h35 mm / approx. 1.7 kg
Protection Category	IP 65



Magnetic Stirrer for Cell Culture, BioMix 1 셀 컬처용 자력 교반기

Cat.No.	Model	Speed Range	Capacity	EA Price ₩
S08-75-075	85001	5 ~ 250 rpm	1 ~ 5000 ml	2,130,200

- Cell Culture나 Tissue Culture용, All Stainless Steel 재질
- 극저속, 고점도, 대용량 및 장시간 사용에서도 안정적인 교반력을 가짐
- Digital Display, Soft Start, Jerk-free Stirring • 2 mag

Cat. No.	S08-75-075
Model	85001
Stirring Point	1
Stirring Volume	5 ~ 5000 ml
Speed Range	5 ~ 250 rpm
Stirring Power	6 W
Power Setting	10-steps (10 ~ 100%)
Material Housing	Stainless Steel
Operation Conditions	-10 up to + 50°C (at 80% Humidity)
Electrical Data	100 ~ 240 V / 50 ~ 60 Hz / 1.5 A
Size / Weight	w180×d238×h38 mm / approx. 1.4 kg
Protection Category	IP 64



4-Point Magnetic Stirrer, MIX 4 MS 4구 자력 교반기, Individually Control

Cat.No.	Model	Speed Range	Capacity	EA Price ₩
S08-75-080	30404	100 ~ 2000 rpm	1 ~ 3000 ml	3,853,100

- 4 Point를 독립적으로 동시에 속도 조절 할 수 있음
- Ultra Flat, Wear-free 자력 교반기, All Stainless Steel 재질
- Digital Display, Soft Start, Jerk-free Stirring
- 2 mag

Cat. No.	S08-75-080
Model	30404
Stirring Point / Distance	4 Point / 140 mm
Stirring Volume	1 ~ 3000 ml
Speed Range	100 ~ 2000 rpm
Stirring Power	28 W (4×7 W)
Material Housing	Stainless Steel
Operation Conditions	-10 up to + 50°C (at 80% humidity)
Electrical Data	100 ~ 240 V / 50 ~ 60 Hz / 1.5 A
Size / Weight	w245×d425×h38 mm / approx. 8.9 kg
Protection Category	IP 65

