

Magnetic Stirrer for Viscous Media 고 점도용 자력 교반기

Cat.No.	Model	Type	Capacity	EA Price ₩
S08-75-481	30001	MIX 1	1 ~ 10,000 ml	1,341,300
S08-75-485	30201	MIX 1 XL	1,000 ml ~ 40,000 ml	2,491,900
S08-75-489	30501	Max MIX	1 L ~ 150 L	5,440,800

- 고점도 및 대용량 시료의 자력 교반에 알맞음
- Digital Display : Speed, Power, Softstart and Interval Setting
- Quick Start and Quick Stop Function
- All Stainless Steel 재질
- 2 mag



Cat. No.	S08-75-481	S08-75-485	S08-75-489
Model	30001	30201	30501
Stirring Point	1		
Stirring Volume	1 ~ 10,000 ml	1,000 ~ 40,000 ml	1 L ~ 150 L
Speed Range	100 ~ 2000 rpm	70 ~ 2000 rpm	100 ~ 1200 rpm
Stirring Power	5 / 10 / 15 / 20 W (4 Step)		40 W
Material Housing	Stainless Steel		
Operation Conditions	-10 °C up to + 50 °C (at 80% Humidity)		
Interval	-	0 ~ 240 sec	-
SoftStart	-	-	Variable, 20 ~ 99 sec
Electrical Data	100 ~ 240 V / 50 ~ 60 Hz / 1.5 A		
Size	w188 × d238 × h38 mm	w225 × d280 × h35 mm	w265 × d320 × h68 mm
Weight	approx. 4.0 kg		approx. 9.3 kg
Protection Category	IP 65	IP 64	

Multiple Magnetic Stirrer, MIX6, MIX15 & MIX15 eco 멀티 자력 교반기

- Ultra Flat, Wear-free 자력 교반기, All Stainless Steel 재질
- Digital Display, Soft Start, Jerk-free Stirring
- 2 mag

6구 멀티 자력 교반기

Cat.No.	Model	Type	Position	Capacity	EA Price ₩
S08-75-192	30006	MIX 6	6	1 ~ 3000 ml	3,348,300

15구 멀티 자력 교반기

Cat.No.	Model	Type	Position	Capacity	EA Price ₩
S08-75-196	30015	MIX 15	15	1 ~ 3000 ml	4,197,000
S08-75-200	30115	MIX 15 eco		1 ~ 1500 ml	2,215,800



Cat. No.	S08-75-192	S08-75-196	S08-75-200
Model	30006	30015	30115
Stirring Point	6	15	
Stirring Point Distance	130 mm	65 mm	
Stirring Volume	1 ~ 3000 ml		1 ~ 1500 ml
Speed Range	100 ~ 2000 rpm		120 ~ 1200 rpm
Stirring Power	5 / 10 / 15 / 20 W (4 Step)		10 W
Material Housing	Stainless Steel		
Operation Conditions	-10 up to + 50°C (at 80% Humidity)		
Electrical Data	100 ~ 240 V / 50 ~ 60 Hz / 1.5 A		
Size	w245 × d425 × h38 mm		
Weight	approx. 8.8 kg		
Protection Category	IP 65		