

PP Handle Beaker PP 핸들 비이커

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- LabPlasti®

Mould 눈금

Cat. No.	Model	용량	Size	EA Price ₩
B05-20-180	15702-0100E	100 ml	Ø61×h72 mm	1,800
B05-20-188	15702-0250E	250 ml	Ø80×h100 mm	1,880
B05-20-204	15702-0500E	500 ml	Ø102×h120 mm	2,250
B05-20-220	15702-1LE	1 L	Ø122×h150 mm	3,010
B05-20-228	15702-2LE	2 L	Ø151×h187 mm	6,010
B05-20-236	15702-3LE	3 L	Ø172×h215 mm	8,650
B05-20-244	15702-5LE	5 L	Ø199×h253 mm	12,020

Blue Printed 눈금

Cat. No.	Model	용량	Size	EA Price ₩
B05-20-280	15702-0100B	100 ml	Ø61×h72 mm	2,000
B05-20-288	15702-0250B	250 ml	Ø80×h100 mm	3,000
B05-20-304	15702-0500B	500 ml	Ø102×h120 mm	3,500
B05-20-320	15702-1LB	1 L	Ø122×h150 mm	4,800
B05-20-328	15702-2LB	2 L	Ø151×h187 mm	9,200
B05-20-336	15702-3LB	3 L	Ø172×h215 mm	12,700
B05-20-344	15702-5LB	5 L	Ø199×h253 mm	18,800



PP Beaker with Handle PP 플라스틱 핸들 비이커

Cat. No.	Model	용량	Size	EA Price ₩	PK Price ₩
B05-46-750	T431040	100 ml	Ø63×h74 mm	2,200	12ea 24,600
B05-46-753	T431050	250 ml	Ø84×h94 mm	3,100	12ea 34,600
B05-46-756	T431060	600 ml	Ø109×h130 mm	4,600	12ea 51,000
B05-46-759	T431070	1 L	Ø129×h150 mm	8,200	6ea 46,000
B05-46-762	T431080	2 L	Ø155×h190 mm	15,000	4ea 56,000
B05-46-765	T431010	3 L	Ø175×h215 mm	15,800	4ea 60,100
B05-46-768	T431090	5 L	Ø205×h265 mm	26,800	2ea 50,000

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Black Graduation
- Tarsons



PMP Beaker with Handle PMP 투명 플라스틱 핸들 비이커

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩	PK Price ₩
B05-46-551	T441040	100 ml	63×74 mm	2,300	12ea 26,100
B05-46-553	T441050	250 ml	84×94 mm	3,800	12ea 44,000
B05-46-555	T441060	600 ml	109×130 mm	6,400	12ea 72,800
B05-46-557	T441070	1 L	129×150 mm	12,000	6ea 68,700
B05-46-559	T441080	2 L	155×190 mm	19,300	4ea 73,700
B05-46-560	T441071	3 L	175×215 mm	29,800	4ea 113,100
B05-46-561	T441090	5 L	205×265 mm	34,200	2ea 65,100

- 재질 : PMP, 0 ~ 150°C 내열, Autoclavable
- 투명 PMP재질로 내열, 내화학성이 뛰어남
- Black Graduation
- Tarsons

