

PFA Beaker with Handle PFA 테프론 반투명 핸들 비이커

Cat. No.	Model	용량	Size	EA Price ₩
B05-24-519	11110	100 ml	∅56.8×h71 mm	64,200
B05-24-523	11111	200 ml	∅69×h92.8 mm	86,100
B05-24-531	11112	500 ml	∅90×h121 mm	176,600
B05-24-535	11113	1 L	∅110×h152 mm	390,800
B05-24-539	11114	2 L	∅136×h198 mm	626,500

- 재질 : PFA, -200 ~ 250°C 내열
- 내부 바닥이 평평하여 자석교반기를 이용한 교반이 원활함
- Universal



PC Handle Beaker 투명 PC 핸들 비이커

Cat. No.	Model	용량	Size	EA Price ₩
B05-517-025	JPC0250P	250 ml	∅70×h100 mm	3,100
B05-517-050	JPC0500P	500 ml	∅100×h125 mm	5,000
B05-517-100	JPC1000P	1 L	∅120×h160 mm	7,900
B05-517-897	JPC2000P	2 L	∅145×h220 mm	11,900
B05-517-900	JPC3000P	4 L	∅180×h265 mm	29,600

- 재질 : PC, -30 ~ 100°C 내열
- Glass와 같은 투명성을 가짐
- LabPlasti®



PE Handle Beaker, 10L 10 L PE 핸들 비이커

Cat. No.	Model	용량	눈금	Size (상/하 ∅ × h)	EA Price ₩
B05-173-201	301204	10 L	500 ml	325/300×266 mm	48,200

- 재질 : HDPE, -40 ~ 85°C 내열, Non-autoclavable
- 손잡이가 2개 부착되어 있어 작업이 용이함
- Universal



Porcelain Beaker 자제 비이커

Cat. No.	Model	용량	Size	EA Price ₩
B05-128-104	227/1	275 ml	∅70×h85 mm	13,500
B05-128-108	227/2	400 ml	∅80×h99 mm	16,300
B05-128-112	227/3	620 ml	∅90×h116 mm	20,500
B05-128-116	227/4	860 ml	∅100×h129 mm	26,100
B05-128-120	227/5	950 ml	∅123×h100 mm	34,200

- 재질 : Porcelain, 1000°C 내열
- Jipo / LabCermic®



Porcelain Beaker with Handle 손잡이형 자제 비이커

Cat. No.	Model	용량	Size	EA Price ₩
B05-128-701	241a/1	100 ml	∅71×h72 mm	22,000
B05-128-703	241a/2	250 ml	∅82×h96 mm	27,300
B05-128-705	241a/3	500 ml	∅97×h131 mm	34,200
B05-128-707	241a/4	1 L	∅127×h160 mm	60,600

- 재질 : Porcelain, 1000°C 내열
- 내부에 눈금이 있어 용량 확인이 가능함
- Jipo / LabCermic®

