

## IKA Stirring Element for Overhead Stirrer 오버헤드 스테리용 임펠라, Stainless Steel AISI 316L

• IKA

### Propeller Stirrer

#### 4 Blade, ①

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-512	0741000	R1342	Ø50 mm	Ø8×L350 mm	95,000
S09-05-518	0741300	R1345	Ø100 mm	Ø8×L540 mm	240,000
S09-05-519	0741300-1	R1345-1	Ø100 mm	Ø8×L350 mm	451,000
S09-05-530	0739000	R2302	Ø150 mm	Ø13×L800 mm	710,000

#### 3 Blade, ②

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-536	1296000	R1381	Ø45 mm	Ø8×L350 mm	610,000
S09-05-542	1295900	R1382	Ø55 mm	Ø8×L350 mm	640,000
S09-05-548	0477700	R1385	Ø140 mm	Ø10×L550 mm	460,000
S09-05-554	0477800	R1388	Ø140 mm	Ø10×L800 mm	570,000
S09-05-556	2343600	R1389	Ø75 mm	Ø8×L350 mm	1,060,000

• R1389는 PTFE Blade & PTFE Coating Shaft임

### Turbine Stirrer, ③

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-566	2332900	R1311	Ø30 mm	Ø8×L350 mm	420,000
S09-05-572	2333000	R1312	Ø50 mm	Ø8×L350 mm	540,000
S09-05-578	2333100	R1313	Ø70 mm	Ø10×L400 mm	640,000

### Dissolver Stirrer, ④

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-584	0513500	R1300	Ø80 mm	Ø8×L350 mm	520,000
S09-05-590	2387900	R1302	Ø100 mm	Ø10×L350 mm	1,110,000
S09-05-596	2746700	R1303	Ø40 mm	Ø8×L350 mm	820,000

### Centrifugal Stirrer, ⑤

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-602	0756900	R1352	Ø60 / 15 mm	Ø8×L350 mm	350,000
S09-05-608	1132700	R1355	Ø100 / 24 mm	Ø8×L550 mm	650,000

### Paddle Stirrer, ⑥

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-620	0757700	R1375	Ø70 mm	Ø8×L550 mm	370,000
S09-05-626	0758000	R1376	Ø150 mm	Ø10×L550 mm	740,000
S09-05-632	0739500	R2311	Ø150 mm	Ø13×L800 mm	1,220,000

### Anchor Stirrer, ⑦

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-638	2022300	R1330	Ø45 mm	Ø8×L350 mm	610,000
S09-05-644	2022400	R1331	Ø90 mm	Ø8×L350 mm	690,000
S09-05-656	2747400	R1333	Ø150 mm	Ø10×L550 mm	840,000

### Moebius Stirrer, ⑧

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-700	20001192	R3000.1	Ø100 mm	Ø8×L565 mm	850,000
S09-05-706	20001195	R3001.1		Ø10×L575 mm	900,000

### Spiral Stirrer, ⑨

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-712	20001203	R3003	Ø50 mm	Ø8×L350 mm	910,000
S09-05-718	20001204	R3003.1	Ø100 mm	Ø10×L550 mm	1,170,000
S09-05-724	20001205	R3003.2	Ø150 mm	Ø8×L800 mm	1,500,000

### Blade Stirrer, ⑩

Cat. No.	Model	Type	Stirring	Shaft Size	EA Price ₩
S09-05-730	20001206	R3004	Ø30 mm	Ø8×L359 mm	650,000
S09-05-736	20001207	R3004.1	Ø50 mm	Ø10×L565 mm	970,000
S09-05-742	20001208	R3004.2	Ø70 mm	Ø10×L819 mm	1,260,000

