

## Stirrer Rod / Impeller, 304 Stainless Steel

### 교반봉 / 임펠러

- 재질 : 304 Stainless Steel
- LK Labkorea

#### Clover-type Impeller, 3-Blade ①

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-111	SR-C040	Ø40 mm	Ø8×L300 mm	21,300
S09-03-118	SR-C060	Ø60 mm	Ø8×L500 mm	24,600
S09-03-125	SR-C080	Ø80 mm	Ø8×L500 mm	29,600
S09-03-132	SR-C100	Ø100 mm	Ø8×L650 mm	31,000

#### Propeller-type Impeller ②

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-140	SR-S070	Ø70 mm	Ø8×L500 mm	30,000
S09-03-148	SR-S100	Ø90 mm	Ø8×L500 mm	82,300
S09-03-156	SR-S120	Ø120 mm	Ø8×L650 mm	86,100

#### Diamond-type Impeller ③

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-170	SR-M090	Ø90 mm	Ø8×L500 mm	18,100

#### Centrifugal-type Impeller

##### Type I ④

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-177	SR-F075	Ø100 mm	Ø8×L500 mm	54,200

##### Type II ⑤

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-179	SR-F080	Ø80 mm	Ø8×L500 mm	60,300
S09-03-185	SR-F100	Ø100 mm	Ø8×L650 mm	104,000

#### Anchor-type Impeller ⑥

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-201	SR-A050	Ø50×h50 mm	Ø8×L500 mm	34,500
S09-03-207	SR-A070	Ø70×h70 mm		35,400
S09-03-209	SR-A080	Ø80×h80 mm		78,100
S09-03-213	SR-A090	Ø90×h90 mm	Ø8×L650 mm	80,000
S09-03-214	SR-A100	Ø100×h100 mm		93,800

#### Round Anchor-type Impeller ⑦

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-222	SR-A150	Ø70 ×h30 mm	Ø8×L500 mm	54,900
S09-03-229	SR-A170	Ø90 ×h40 mm		62,700
S09-03-236	SR-A190	Ø110 ×h50 mm	Ø8×L650 mm	78,400

#### 반달형 Impeller

##### 고정형 ⑧

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-274	SR-H1080	Ø80 ×h30 mm	Ø8×L500 mm	50,200
S09-03-281	SR-H1100	Ø100 ×h40 mm		54,900
S09-03-288	SR-H1120	Ø120 ×h50 mm	Ø8×L650 mm	62,700

- 하부 Blade가 Impeller에 고정되어 있음

##### 회전형 ⑨

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-653	SR-TH060	Ø60 ×h19 mm	Ø8×L500 mm	59,800
S09-03-660	SR-TH080	Ø80 ×h19 mm		67,600
S09-03-667	SR-TH100	Ø100 ×h23 mm		67,600

- PTFE Blade & PTFE Coating Shaft임

#### Turbin-type Impeller ⑩

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-315	SR-T040	Ø40 mm	Ø8×L300 mm	58,500
S09-03-322	SR-T070	Ø70 mm	Ø8×L500 mm	65,000

