


Stirring Rod-Impeller 에서 계속...

Centrifugal-type Impeller


Type I ⑤

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-579	SR-TF080	Ø80 mm	Ø8×L500 mm	34,800
S09-03-585	SR-TF100	Ø120 mm	Ø8×L650 mm	55,700



Type II ⑥

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-450-302	005.430.8.2	Ø40 mm	Ø8×L300 mm	102,800
S09-450-310	005.440.8.2		Ø8×L400 mm	103,600
S09-450-318	005.450.8.2	Ø50 mm	Ø8×L500 mm	105,100
S09-450-326	005.455.10.2	Ø70 mm	Ø10×L550 mm	124,000
S09-450-334	005.0465.10.2		Ø10×L650 mm	125,100
S09-450-342	005.0475.10.2		Ø10×L750 mm	126,400
S09-450-350	005.0475.12.2	Ø80 mm	Ø12×L750 mm	163,700



• Cowie

Type III ⑦

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-238-265	18900078	Ø85 mm	Ø8×L350 mm	61,900



반달형 Impeller

고정형 ⑧


Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-140-111	SR-T1080	Ø80 ×h30 mm	Ø8×L500 mm	115,200
S09-140-115	SR-T1100	Ø100 ×h40 mm		112,500
S09-140-119	SR-T1120	Ø120 ×h50 mm	Ø8×L650 mm	143,700



• 하부 Blade가 Impeller에 고정되어 있음

회전형 ⑨


Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-253	SR-H060	Ø60 ×h19 mm	Ø8×L500 mm	32,500
S09-03-260	SR-H080	Ø80 ×h19 mm		45,500
S09-03-267	SR-H100	Ø100 ×h23 mm		32,500



• Blade가 360° 회전 가능하도록 Rod 중간에 홈이 있음

Anchor-type Impeller ⑩


Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-601	SR-TA060	Ø60×h60 mm	Ø8×L500 mm	72,800
S09-03-607	SR-TA080	Ø80×h60 mm		81,600
S09-03-613	SR-TA100	Ø100×h100 mm	Ø8×L650 mm	100,000



Round Anchor-type Impeller

Type I ⑪

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-622	SR-TA160	Ø60 ×h30 mm	Ø8×L500 mm	59,800
S09-03-629	SR-TA180	Ø80 ×h40 mm		67,600
S09-03-636	SR-TA200	Ø100 ×h50 mm	Ø8×L650 mm	107,900



Type II ⑫

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-450-703	005.330.2	Ø80 mm	Ø6×L300 mm	71,600
S09-450-710	005.340.2		Ø6×L400 mm	72,400
S09-450-717	005.350.2		Ø6×L500 mm	73,400
S09-450-719	005.330.82		Ø8×L300 mm	89,800
S09-450-721	005.340.82		Ø8×L400 mm	90,500
S09-450-723	005.350.82		Ø8×L500 mm	91,700
S09-450-724	005.355.10.2	Ø100 mm	Ø10×L550 mm	118,800
S09-450-731	005.0365.10.2	Ø140 mm	Ø10×L650 mm	120,000
S09-450-738	005.0375.10.2		Ø10×L750 mm	121,300



• Cowie • 005.0XXX.10.2는 L50 × Ø6.5mm Stainless Steel Core end

Turbin-type Impeller ⑬

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-03-715	SR-TF040	Ø40 mm	Ø8×L300 mm	153,400
S09-03-722	SR-TF070	Ø70 mm	Ø8×L500 mm	171,600

