

"BOLA" Stirring Element for Overhead Stirrer 오버헤드 스티어링 임펠러

- 재질 : PTFE, -200 ~ 260°C 내열
- BOLA

Gassing Stirrer, ①

Cat. No.	Model	Shaft dia x Length	Chucking dia	Dimension(AxBxCxD)	EA Price ₩
S09-460-421	C488-08	∅ 10 X 484 mm	∅ 8.0 mm	72 x 12 x 20 x 187 mm	935,900
S09-460-425	C488-14	∅ 10 X 559 mm		72 x 12 x 20 x 272 mm	931,100
S09-460-429	C488-20	∅ 10 X 657 mm		72 x 12 x 20 x 387 mm	980,600

Stirrer Shaft, Fan-Shaped, ②

Cat. No.	Model	Shaft dia x Length	Chucking dia	Joint	Dimension(AxB)	EA Price ₩
S09-460-391	C382-02	∅ 8 x 300 mm	∅ 6.5 mm	29/32	24 x 35 mm	556,100
S09-460-395	C382-06			45/40	38 x 45 mm	828,300
S09-460-399	C382-08	∅ 8 x 450 mm	∅ 8.0 mm	60/46	53 x 55 mm	726,500
S09-460-407	C382-14	∅ 10 x 600 mm				847,500

Double Impulse Stirrer Shaft, ③

Cat. No.	Model	Shaft dia x Length	Chucking dia	Blade	Dimension(AxBxC)	EA Price ₩
S09-460-251	C391-18	∅ 10 x 600 mm	∅ 8.0 mm	150mm	140 x 34 x 19 mm	1,487,100
S09-460-255	C391-28	∅ 16 x 800 mm	∅ 14.0 mm		140 x 34 x 19 mm	1,880,700

Slip-on Baffle, ④




Cat. No.	Model	Joint	Width	Shaft dia x Fitting	EA Price ₩
S09-460-371	C490-10	19/26	15 mm	∅ 8 x L125 mm	63,600
S09-460-375	C490-12	29/32	23 mm		41,800

Micro Surface Stirrer Shaft, ⑤

Cat. No.	Model	Shaft dia x Length	Chucking dia	Dimension(AxB)	EA Price ₩
S09-460-151	C486-08	∅ 3.5 x 120 mm	∅ 2.5 mm	8 x 8 mm	235,700
S09-460-155	C486-12	∅ 3.5 x 180 mm		12 x 12 mm	251,800
S09-460-159	C486-16	∅ 4.0 x 200 mm	∅ 3.0 mm	14 x 14 mm	261,300
S09-460-163	C486-20			16 x 16 mm	250,200

Mini-Propeller Stirrer Shaft, ⑥

Cat. No.	Model	Shaft dia x Length	사용 Suitable	Dimension(AxBxC)	EA Price ₩
S09-460-181	C482-08	∅ 6 x 200 mm	NS 29/32	25 x 8 x 2mm	212,800
S09-460-185	C482-12	∅ 6 x 350 mm		259,500	
S09-460-189	C482-20	∅ 6 x 200 mm	NS 45/40	40 x 12 x 2mm	265,100
S09-460-193	C482-24	∅ 6 x 350 mm			241,100



- 316 Stainless Steel Impeller는 1040 Page 참조

