

## Adjustable PTFE Blade Only 이동식 테프론 블레이드

- 재질 : PTFE, Glass Filled PTFE Lock Nut
- Cowie

### 4 Blade Angled Type, 45°

Cat. No.	Model	Blade Size	Shaft Size	EA Price ₩
S09-450-555	005.608040	Ø40 mm	Ø8 mm	81,200
S09-450-558	005.608070	Ø70 mm		109,600
S09-450-561	005.608090	Ø90 mm		109,000
S09-450-564	005.618060	Ø60 mm	Ø10 mm	87,300
S09-450-567	005.618090	Ø90 mm		131,100

- Blade의 하부쪽으로 흐름을 만들어 줌

### 4 Blade Flat Type, 135°

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-450-587	005.608041	Ø40 mm	Ø8 mm	102,400
S09-450-591	005.608071	Ø70 mm		113,300
S09-450-594	005.608091	Ø90 mm		102,400

- Blade의 상부쪽으로 흐름을 만들어 줌

### 2 Blade Angled Type, 90°

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-450-571	005.708070	Ø70 mm	Ø8 mm	68,700
S09-450-573	005.710070		Ø10 mm	68,700
S09-450-577	005.710100	Ø100 mm	Ø10 mm	82,700

### Turbine Blade Type

Cat. No.	Model	Blade Size	Rod Size	EA Price ₩
S09-450-611	005.90850	Ø50 mm	Ø8 mm	98,600
S09-450-616	005.90875	Ø75 mm		125,100
S09-450-621	005.91075		Ø100 mm	Ø10 mm
S09-450-626	005.910100	165,900		
S09-450-631	005.91275	Ø75 mm		Ø12 mm
S09-450-636	005.912100	Ø100 mm	Ø12 mm	174,800
S09-450-641	005.916100			183,300
S09-450-646	005.916150	Ø150 mm	Ø16 mm	291,000




4 Blade, 45° Flat



4 Blade, 135° Angle



2 Blade, 90° Angle



Turbine Blade Type

## PTFE Stirring Rod Only PTFE 교반봉

- 재질 : PTFE, -200 ~ 280°C 내열
- 불투명 PTFE 테프론 재질로 최상의 내화학성 가짐
- Steel Core Type can be bent into a permanent shape
- Cowie

### Steel Core Type ①

Cat. No.	Model	Size	EA Price ₩
S09-450-401	003.100	Ø6×L100 mm	8,800
S09-450-403	003.150	Ø6×L150 mm	9,600
S09-450-405	003.200	Ø6×L200 mm	10,600
S09-450-407	003.250	Ø6×L250 mm	11,800
S09-450-409	003.300	Ø6×L300 mm	14,100
S09-450-411	003.350	Ø6×L350 mm	17,900
S09-450-413	003.400	Ø6×L400 mm	20,200

### Solid PTFE Type, Tapered End ②

Cat. No.	Model	Size	EA Price ₩
S09-450-415	003.1100	Ø8×L100 mm	5,100
S09-450-417	003.1150	Ø8×L150 mm	6,400
S09-450-419	003.1200	Ø8×L200 mm	8,000
S09-450-421	003.1250	Ø8×L250 mm	9,400
S09-450-423	003.1300	Ø8×L300 mm	10,700

