

## Teflon Halfmoon shape Impeller 반달형 교반봉, Blade 교체식



- 재질 : PTFE, -200 ~ 280°C 내열
- Cowie

### PTFE Rod Only ①

Cat. No.	Model	Type	Size	EA Price ₩
S09-450-127	005.540.8.1	Rod only	Ø8×L400 mm	64,000
S09-450-129	005.550.8.1		Ø8×L500 mm	65,200

### PTFE Blade Only

Cat. No.	Model	Type	Size	EA Price ₩
S09-450-201	002.052.1	Square End ②	w52×h14 mm	4,300
S09-450-205	002.076.1		w76×h19 mm	5,600
S09-450-209	002.090.1		w90×h28 mm	7,700
S09-450-213	002.1065.1	Plain End ③	w65×h25 mm	5,900
S09-450-217	002.1075.1		w75×h25 mm	6,300
S09-450-221	002.1105.1		w105×h25 mm	7,100
S09-450-225	002.1125.1		w125×h25 mm	7,400
S09-450-229	002.1150.1		w150×h25 mm	8,200



## Glass Stirrer Bearing Set 유리 교반봉 세트, LukeGL®

- LK Labkorea

Cat. No.	Model	Joint	Blade (w×h)	Shaft (od×L)	EA Price ₩
S09-132-215	2060-01	14/20	33×11mm	6×305 mm	50,600
S09-132-223	2068-01	24/40	48×19 mm	10×445 mm	60,800
S09-132-231	2068-02		60,800		
S09-132-239	2068-03	29/42	75×19 mm	64,100	
S09-132-246	2068-04	34/45		67,500	
S09-132-254	2068-05	45/50	110×24 mm	10×515 mm	84,300

- 구성 : Glass Bearing, Glass Shaft, PTFE Blade

### Glass Stirrer Bearing Only

Cat. No.	Model	Joint	EA Price ₩
S09-132-277	2060-10	14/20	20,000
S09-132-284	2071-01	24/40	32,900
S09-132-291	2071-02	29/42	33,600
S09-132-298	2071-04	34/45	36,500
S09-132-305	2071-05	45/50	33,700

### Glass Stirring Shaft Only

Cat. No.	Model	O.D	Length	EA Price ₩
S09-132-555	2060-11	Ø6 mm	305 mm	26,700
S09-132-563	2075-01	Ø10 mm	445 mm	30,400
S09-132-571	2075-02		515 mm	33,700
S09-132-589	2075-03		590 mm	37,100

### PTFE Blade Only ①

Cat. No.	Model	Size (mm)	사용 Joint	EA Price ₩
S09-132-665	2080-00	w33×h11, 3t	14/20이상	8,500
S09-132-669	2080-01	w48×h19, 3t	24/40이상	9,700
S09-132-676	2080-02	w60×h19, 3t		9,900
S09-132-683	2080-03	w75×h19, 3t	29/42이상	10,200
S09-132-690	2080-04	w75×h23, 3t		10,300
S09-132-695	2080-05	w110×h24, 3t	34/45이상	11,500
S09-132-699	2080-06	w110×h28, 3t		12,200

### Jacketed Glass Stirrer Bearing Only ②

Cat. No.	Model	Joint	사용 Shaft	EA Price ₩
S09-132-312	2072-01	24/40	Ø10 mm	42,200
S09-132-319	2072-02	29/42		47,200
S09-132-326	2072-03	34/45		50,600
S09-132-333	2072-04	45/50		67,500

