

## Glass Stopcock, PTFE-Plug 유리 스탑 콕크, 테프론 플러그

- Borosilicate Glass 3.3
- Hose Nipple 형
- LukeGL®

### 일자형 유리 스탑 콕크 / 2-way Glass Stopcock

Cat. No.	Model	Bore	Stopcock NS	EA Price ₩
S11-11-111	GC-025	2.5 mm	14.5	33,700
S11-11-119	GC-040	4.0 mm	18.8	50,700
S11-11-127	GC-060	6.0 mm	21.5	72,400
S11-11-135	GC-080	8.0 mm	24	121,400

### 삼방 유리 스탑 콕크 / 3-way Glass Stopcock


Cat. No.	Model	Bore	Stopcock NS	EA Price ₩
S11-11-211	GC-T025	2.5 mm	14.5	89,500
S11-11-219	GC-T040	4.0 mm	18.8	106,400
S11-11-227	GC-T060	6.0 mm	29.2	177,300



## PP Stopcock, PTFE-Plug 플라스틱 스탑 콕크, 테프론 플러그

Cat.No.	Model	Type	Bore	EA Price ₩
S11-53-265	890010	일자형	3.5 mm	22,900
S11-53-284	890020	삼방		27,700

- 재질 : PP(Body) / PTFE(Plug)
- Hose Connection : Ø7 ~ 10 mm
- Tarsons



## Plastic Stopcock 플라스틱 스탑 콕크, 테프론 플러그

- 재질 : PP / PE, 5 ~ 40°C Operating Temp. PVDF, 0 ~ 40°C Operating Temp
- Max. Pressure : 1 bar
- Wenk Labtec

### 일자형 플라스틱 스탑 콕크 / 2-way Plastic Stopcock

Cat. No.	Model	Bore	재질	Tube id	EA Price ₩
S11-20-011	SWQ212	4 mm	PVDF	5 ~ 7 mm	25,500
S11-20-017	SWQ214	6 mm		7 ~ 9 mm	30,200
S11-20-023	SWQ216	8 mm		9 ~ 11 mm	34,600
S11-20-025	SWQ217	8 mm	PE	11 ~ 13 mm	34,600
S11-20-311	SWQ200	4 mm		5 ~ 7 mm	15,300
S11-20-317	SWQ202	6 mm		7 ~ 9 mm	17,600
S11-20-323	SWQ204	8 mm		9 ~ 11 mm	20,100
S11-20-325	SWQ205	8 mm		11 ~ 13 mm	20,100



### 삼방 플라스틱 스탑 콕크 / 3-way Plastic Stopcock

Cat. No.	Model	Bore	재질	Tube id	EA Price ₩
S11-20-111	SWQ218	4 mm	PVDF	5 ~ 7 mm	39,800
S11-20-117	SWQ220	6 mm		7 ~ 9 mm	44,900
S11-20-123	SWQ224	8 mm		9 ~ 11 mm	48,400
S11-20-125	SWQ225	8 mm		11 ~ 13 mm	48,400
S11-20-411	SWQ206	4 mm	PE	5 ~ 7 mm	25,400
S11-20-417	SWQ208	6 mm		7 ~ 9 mm	30,000
S11-20-423	SWQ210	8 mm		9 ~ 11 mm	34,500
S11-20-425	SWQ211	8 mm		11 ~ 13 mm	34,500

- 윈터치 Ball Valve는 291 Page 참조
- PFA Needle Valve는 293 Page 참조
- Stainless Steel Needle Valve는 299 Page 참조
- 초자가공용 Glass & PTFE Cock는 551 ~ 552 Page 참조
- 고진공용 Stopcock는 552 ~ 553 Page 참조
- J-Young® 테프론 Needle Valve는 552 ~ 553 Page 참조

