

PFA Stopcock, Screw-type

테프론 스탑콕크


- 재질 : PFA (Body), PTFE (Plug)
- Temp. Range : -18 ~ 121°C
- Pressure Rating : 50 psi up to 60°C
- Bore : Ø3.2 mm
- Tommyheco

일자형 PFA 테프론 스탑 콕크 / 2-way PFA Stopcock

Cat. No.	Model	Connection	EA Price ₩
S11-255-555	PVT2-T2	od1/8" Tubing	188,200
S11-255-562	PVT4-T4	od1/4" Tubing	216,000
S11-255-569	PVT4-2	od1/4" Tubing, 1/8" MNPT	165,800
S11-255-576	PVT4-4	od1/4" Tubing, 1/4" MNPT	180,700

삼방 PFA 테프론 스탑 콕크 / 3-way PFA Stopcock

Cat. No.	Model	Connection	EA Price ₩
S11-34-024	PVT4-T4-3	od1/4" Tubing	286,900



PTFE Stopcock

테프론 스탑콕크


- 재질 : PTFE (Body), PP (Handle)
- 사용조건 : 5 mmHg (진공), 1 bar (압력)
- Cowie

일자형 테프론 스탑 콕크 / 2-way PTFE Stopcock

Cat. No.	Model	Bore	Connection	EA Price ₩
S11-450-233	016.702.5.2	2 mm	od 4.5 mm	37,600
S11-450-237	016.702.2		od 6 mm	37,600
S11-450-240	016.703.7.2	3 mm	od 6.8 mm	39,500
S11-450-245	016.703.2	3 mm	od 8 mm	39,500
S11-450-247	016.704.9.2	4 mm	od 9 mm	39,800
S11-450-253	016.704.2	4 mm	od 10 mm	40,800
S11-450-257	016.705.2	5 mm	od 11 mm	65,200

삼방 테프론 스탑 콕크 / 3-way PTFE Stopcock

Cat. No.	Model	Bore	Connection	EA Price ₩
S11-450-273	016.802.5.2	2 mm	od 4.5 mm	43,400
S11-450-277	016.802.2		od 6 mm	43,700
S11-450-280	016.803.7.2	3 mm	od 6.8 mm	45,900
S11-450-285	016.803.2		od 8 mm	47,200
S11-450-287	016.804.9.2	4 mm	od 9 mm	48,100
S11-450-293	016.804.2		od 10 mm	48,500
S11-450-295	016.805.2	5 mm	od 11 mm	62,300



PTFE Stopcock, Screw-type

테프론 스탑콕크

- 재질 : PTFE (Body), PP (Handle)
- 사용조건 : 5 mmHg (진공), 1 bar (압력)
- Cowie

일자형 테프론 스탑 콕크 / 2-way PTFE Stopcock

Cat.No.	Model	Bore	Connection	EA Price ₩
S11-450-261	016.1202.6.2	2 mm	od6 mm	44,500
S11-450-269	016.1204.2	4 mm	od8 mm	51,100

삼방 테프론 스탑 콕크 / 3-way PTFE Stopcock

Cat.No.	Model	Bore	Connection	EA Price ₩
S11-450-301	016.1302.6.2	2 mm	od6 mm	61,600
S11-450-309	016.1304.2	4 mm	od8 mm	67,000

