

## PTFE Stopcock with Hose Connection

### PTFE 테프론 호스용 스탑콕크

- 재질 : PTFE, PP, 0 °C to 110°C 내열
- 사용 조건 : Max. 2 bar
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE만 접촉 되어짐
- BOLA

#### 2-Way Stopcock, 일자형

Cat. No.	Model	hose connectors	Tubing id	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-011	E650-03	od 4.5 mm	4 mm	Ø 1.5 mm	60 x 22 x 53 mm	176,300
S11-460-015	E650-06	od 6.8 mm	6 mm	Ø 3.0 mm		126,600
S11-460-019	E650-09	od 9.0 mm	8 mm	Ø 4.0 mm		167,700
S11-460-023	E650-12	od 11.0 mm	10 mm	Ø 6.0 mm	85 x 35 x 69 mm	192,500



2-Way

#### 3-Way Stopcock, 일자형

Cat. No.	Model	hose connectors	Tubing id	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-051	E650-50	od 4.5 mm	4 mm	Ø 1.0 mm	60 x 41 x 53 mm	304,600
S11-460-055	E650-53	od 6.8 mm	6 mm	Ø 2.0 mm		278,700
S11-460-059	E650-56	od 9.0 mm	8 mm	Ø 3.0 mm		295,700
S11-460-063	E650-59	od 11.0 mm	10 mm	Ø 4.0 mm	85 x 60 x 69 mm	348,200



3-Way

## PTFE Stopcock

### PTFE 테프론 스탑콕크

- 재질 : PTFE, 0 °C to 110°C 내열
- 사용 조건 : Max, 2 bar
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE만 접촉 되어짐
- BOLA

#### 2-Way Stopcock, 일자형

Cat. No.	Model	Thread	Tubing od	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-101	E652-02	GL 14	4 mm	Ø 2 mm	59 x 22 x 53mm	142,500
S11-460-105	E652-04		6 mm			145,800
S11-460-109	E652-06	GL 18	8 mm	Ø 5 mm	74 x 35 x 69mm	210,100
S11-460-113	E652-08		10 mm			211,000



2-Way

#### 3-Way Stopcock, T-Bore Shape

Cat. No.	Model	Thread	Tubing od	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-131	E654-02	GL 14	4 mm	Ø 1.5 mm	59 x 41 x 53mm	187,600
S11-460-135	E654-04		6 mm			185,700
S11-460-139	E654-06	GL 18	8 mm	Ø 3.5 mm	74 x 54 x 69mm	283,000
S11-460-143	E654-08		10 mm			285,100



3-Way

## GL Stopcock

### PTFE 테프론 GL 스크류 스탑콕크

- 재질 : PTFE, 0 °C to 110 °C 내열
- 사용 조건 : Max. 6 bar (Pressure) 또는 Vacuum
- 시료는 단지 PTFE 재질에만 접촉 되어짐
- BOLA

#### 2-Way 일자형


Cat. No.	Model	Thread	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-201	E684-14	GL 14	Ø 4 mm	54 x 20 x 38mm	139,000
S11-460-205	E684-18	GL 18	Ø 6 mm	64 x 30 x 45mm	153,600
S11-460-209	E684-25	GL 25	Ø 8 mm	78 x 40 x 57mm	289,500



2-Way

#### 3-Way, L-Bore Shape

Cat. No.	Model	Thread	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-231	E686-14	GL 14	Ø 4 mm	64 x 47 x 43mm	220,000
S11-460-235	E686-18	GL 18	Ø 6 mm	74 x 57 x 57mm	299,800
S11-460-239	E686-25	GL 25	Ø 6 mm	78 x 59 x 57mm	440,000



3-Way, L-Bore Shape

3-Way, T-Bore Shape

#### 3-Way, T-Bore Shape

Cat. No.	Model	Thread	Bore Dia	외부 Size (L x D x H mm)	EA Price ₩
S11-460-251	E688-14	GL 14	Ø 4 mm	74 x 57 x 57mm	381,300
S11-460-255	E688-18	GL 18	Ø 4 mm	74 x 57 x 57mm	315,800
S11-460-259	E688-25	GL 25	Ø 6 mm	88 x 69 x 57mm	432,000