

PTFE Teflon Joint Stopper 기본형 테프론 죠인트 마개

- Extraction Nut를 가지고 있어 쉽게 제거
- Cowie

Medium Length Type

Cat. No..	Model	Joint Size	EA Price ₩
S12-450-355	010.210	10/19	12,800
S12-450-362	010.214	14/23	13,300
S12-450-369	010.219	19/26	16,300
S12-450-376	010.224	24/29	21,700
S12-450-383	010.229	29/32	24,800
S12-450-390	010.234	34/35	33,300

Full Length Type

Cat. No..	Model	Joint Size	EA Price ₩
S12-450-455	010.110	10/30	14,400
S12-450-462	010.114	14/35	20,400
S12-450-469	010.119	19/38	24,400
S12-450-476	010.124	24/40	25,700
S12-450-483	010.129	29/42	33,100
S12-450-490	010.134	34/45	34,900



PTFE Teflon Tight Joint Stopper 테프론 죠인트 마개, Self Releasing

- Cowie

Medium Length Type

Cat. No..	Model	Joint Size	EA Price ₩
S12-450-069	010.1219	19/26	26,000
S12-450-076	010.1224	24/29	34,900
S12-450-083	010.1229	29/32	42,400
S12-450-090	010.1234	34/35	44,700
S12-450-095	010.1245	45/40	79,000

Full Length Type

Cat. No..	Model	Joint Size	EA Price ₩
S12-450-169	010.1119	19/38	29,600
S12-450-176	010.1124	24/40	38,700
S12-450-190	010.1134	34/45	48,900
S12-450-195	010.1145	45/50	86,400



PP Test Tube Stopper, Kimble® 플라스틱 튜브 마개

Cat. No..	Model	사용 Tube od	색상	PK Price ₩	CS Price ₩
S12-55-111	73660-13	Ø13 mm	Natural	100ea 19,900	10pk 186,000
S12-55-114	73662-13		Red	100ea 25,700	10pk 239,400
S12-55-117	73660-16	Ø16 mm	Natural	100ea 21,300	10pk 198,600
S12-55-120	73662-16		Red	100ea 26,600	10pk 248,200
S12-55-123	73660-18	Ø18 mm	Natural	100ea 23,700	10pk 221,100
S12-55-126	73662-18		Red	100ea 33,500	10pk 312,400
S12-55-129	73660-20	Ø20 mm	Natural	100ea 40,200	5pk 187,600
S12-55-132	73662-20		Red	100ea 52,700	5pk 246,000
S12-55-135	73660-25	Ø25 mm	Natural	100ea 62,000	5pk 289,600
S12-55-138	73662-25		Red	100ea 72,300	5pk 337,800
S12-55-141	73660-38	Ø38 mm	Natural	125ea 193,000	2pk 360,700

- 재질 : PP, 0 ~ 135°C 내열, Autoclavable
- 요청에 의해 Yellow, Green, Blue 색상도 공급 가능
- DWK Lifescience / Kimble®



73662-25



73660-38