

Vacuum Sublimation Apparatus 진공 서브리메이션 장치

Cat.No.	Model	용량	Flange	전체 Size (mm)	EA Price ₩
S22-132-321	3038-01	150 ml	80 mm	215×h285	1,086,100

- O-ring Sealed Flange를 가진 장치로 Grease Free
- Dry Ice, Acetone 또는 액체 질소를 이용하여 Condensing
- Hotplate나 Hotplate Stirrer와 같이 사용 가능
- Chemglass

구성 Component

Cat.No.	Model	Description	EA Price ₩
S22-132-532	3038-02	Lower Flask Only	231,200
S22-132-538	3038-03	Top Condenser Only	867,000
S22-132-544	139-03	Clamp, 80 mm Flange	151,900
S22-132-550	305-339	Viton O-ring, #339	18,000
S22-132-556	961-02	0 ~ 8 mm Plug and Control Knob	78,000

Optional Component

Cat.No.	Model	Description	EA Price ₩
S22-132-580	309-339	O-ring, Perfluoro, #339	82,000

- 내화학성이 우수함



Sublimation Apparatus 서브리메이션 장치, LukeGL®

Cat.No.	Model	용량	Joint	EA Price ₩
S22-130-332	SA-M03	3 ml	29/32	116,100
S22-130-340	SA-M12	12 ml	34/45	201,300

- Screw-Thread Joint와 2.5 mm Bore Glass Stopcock
- LukeGL®

구성 Component for SA-M03

Cat.No.	Model	Description	EA Price ₩
S22-130-666	SA-M031	Sublimation Head	67,500
S22-130-678	SA-M032	Sublimation Reservoir	33,700
S22-130-690	SA-M033	Viton O-ring	1,200
S22-130-702	SA-M034	Optional Teflon Coating O-ring	15,200

구성 Component for SA-M12

Cat.No.	Model	Description	EA Price ₩
S22-130-672	SA-M121	Sublimation Head	84,300
S22-130-684	SA-M122	Sublimation Reservoir	50,600
S22-130-696	SA-M123	Viton O-ring	1,600
S22-130-708	SA-M124	Optional Teflon Coating O-ring	15,200



Micro Sublimation Apparatus 마이크로 서브리메이션 장치

Cat.No.	Model	용량	Joint	EA Price ₩
S22-132-357	3030-01	3 ml	24/25	291,900

- Screw-Thread Joint와 2.5 mm Bore Glass Stopcock
- O-ring Seal Joint
- Chemglass

Replacement Part

Cat. No.	Model	Description	EA Price ₩
S22-132-364	3030-02	Lower Flask Only	78,000
S22-132-371	3030-03	Cold Finger Condenser Only	224,900

