

## Disposable PP Syringe 일회용 PP 주사기


- 재질 : PP 몸체 및 PE Plunger (FDA Registerd)
- Plunger에 고무 가스켓이 없는 내화학용 주사기
- 개별 멸균 포장, 주사바늘은 별도 구입 • 유기용매 사용에 적합
- Laboratory & Industrial Use Only
- HSW (Henke Sass Wolf)

### Luer Tip, 기본형

Cat. No.	Model	용량	눈금	PK Price ₩	
S15-07-008	A1	1 ml	0.01 ml	100ea	15,200
S15-07-016	A2	2/3 ml	0.1 ml	100ea	10,300
S15-07-024	A5	5 ml	0.2 ml	100ea	16,300
S15-07-032	A10	10 ml	0.5 ml	100ea	25,200
S15-07-040	A20	20 ml	1 ml	100ea	43,700
S15-07-048	A30	30 ml	1 ml	50ea	54,100
S15-07-056	A50	50 ml	2 ml	30ea	37,700

### Luer-Lock Tip

Cat. No.	Model	용량	눈금	PK Price ₩	
S15-07-064	AL2	2/3 ml	0.1 ml	100ea	27,300
S15-07-072	AL5	5 ml	0.2 ml	100ea	34,800
S15-07-080	AL10	10 ml	0.5 ml	100ea	41,500
S15-07-088	AL20	20 ml	1 ml	100ea	55,600
S15-07-096	AL30	30 ml	1 ml	50ea	57,200
S15-07-104	AL50	50 ml	2 ml	30ea	48,800



## Disposable Syringe 일회용 주사기

Cat. No.	Model	용량	눈금	주사침	PK Price ₩	
S15-27-125	DS001	1 ml	0.01 ml	26G, 2/1"	100ea	8,500
S15-27-132	DS002	2 ml	0.1 ml	23G, 1"	100ea	8,400
S15-27-139	DS003	3 ml	0.1 ml		100ea	8,100
S15-27-146	DS005	5 ml	0.2 ml	22G, 1 1/4"	100ea	9,800
S15-27-153	DS010	10 ml	0.2 ml		100ea	13,300
S15-27-160	DS020	20 ml	1 ml	21G, 1 1/4"	50ea	11,200
S15-27-167	DS030	30 ml	1 ml		50ea	16,800
S15-27-174	DS050	50 ml	2 ml		25ea	10,800

- 국산 주사기로 주사침이 기본 포함
- Laboratory & Industrial Use Only
- LK Labkorea



## Matched Numbered Glass Syringe 고급형 유리 주사기, Metal Luer-Lock Tip

Cat. No.	Model	용량	눈금	EA Price ₩
S15-205-317	2501	1 ml	0.05 ml	124,200
S15-205-338	2717	5 ml	0.2 ml	117,300
S15-205-345	2527	10 ml		122,000
S15-205-352	2537	20 ml	1 ml	138,600
S15-205-359	2547	30 ml		227,600
S15-205-366	2557	50 ml	2 ml	293,500
S15-205-373	2577	100 ml	4 ml	956,300
S15-205-245	2733	250 ml	10 ml	1,962,000

- 각 주사기의 Barrel과 Plunger는 고유 넘버를 가지고 있음
- 정확한 기록을 요구하는 Critical Application용
- 기본형보다 Higher Pressure에 견딜수 있음
- Cadence Science

