

Gas Tight Syringe for Manual Use 가스 샘플 투입용 주사기, SS Plunger

- 주용도 : Sample Dosage for Gas Chromatography, HPLC Valve, Autosampler
- Plunger : Stainless Steel
- 요청에 의해 SS Plunger Type 에 Gold Coating 이 가능함
- Glass Barrel : Borosilicate 3.3
- Needle Length : 51 mm (A Type)
- **Setonic**

Fixed Needle Type


Cat. No.	Model	용량	Plunger Type	Needle		EA Price ₩
				od (Gauge)	길이	
S14-165-401	3030205	5 ul	SS	0.47 mm (G26S)	51 mm	37,000
S14-165-405	3030206	10 ul				31,700
S14-165-409	3030207	25 ul				38,900
S14-165-413	3030208	50 ul				37,000
S14-165-417	3030209	100 ul				44,000
S14-165-421	3030210	250 ul				47,600

Removable Needle Type

Cat. No.	Model	용량	Plunger Type	Needle		EA Price ₩
				od (Gauge)	길이	
S14-165-451	3030389	5 ul	SS	0.47 mm (G26S)	51 mm	60,000
S14-165-455	3030298	10 ul				51,200
S14-165-459	3030299	25 ul		44,000		
S14-165-463	3030300	50 ul		54,600		
S14-165-467	3030301	100 ul		60,000		
S14-165-471	3030302	250 ul		63,500		
S14-165-475	3030303	500 ul		63,500		
				0.72 mm (G22)		

Replacement Needle

Cat. No.	Model	용량	Needle		PK Price ₩	
			od (Gauge)	길이		
S14-165-901	3050072	5 ul	0.47 mm (G26S)	51 mm	3ea	39,900
S14-165-905	3050022	10 ul			3ea	28,300
S14-165-909	3050179	25 uil ~ 10 ml	0.72 mm (G22)		3ea	28,300
S14-165-913	3050036		0.72 mm (G22S)		3ea	28,300



Fixed Needle Type



Removable Needle Type




Replacement Needle

Gas Tight Syringe for Manual Use 가스 샘플 투입용 주사기, PTFE Plunger

- 액상 및 기체상 시료에 모두 사용
- 테프론 팁이 유리 재질의 몸통과 빈틈없이 밀착하여 Tight한 Sealing 상태를 만들어줌
- 용도 : Sample 및 Aggressive Fluid의 Dosage
- Glass Barrel : Borosilicate 3.3
- Needle Length : 51 mm (A Type)
- **Setonic**

Fixed Needle Type

Cat. No.	Model	용량	Plunger Type	Needle		EA Price ₩
				od (Gauge)	길이	
S14-165-301	3030241	5 ul	PTFE	0.47 mm (G26S)	51 mm	51,200
S14-165-305	3030242	10 ul				42,300
S14-165-309	3030243	25 ul		45,900		
S14-165-313	3030244	50 ul		0.72 mm (G22S)		49,300
S14-165-317	3030245	100 ul		52,900		
S14-165-321	3030246	250 ul		0.72 mm (G22)		56,500
S14-165-325	3030247	500 ul				52,900



Fixed Needle Type

다음 페이지에서 계속...