

Manual HPLC Syringe HPLC 용 주사기

- 주용도 : Sample Dosage for Gas Chromatography, HPLC Valve, Autosampler
- Manual Using Waters UK6 Valve
- Glass Barrel : Borosilicate 3.3
- Needle Length : 50 mm (B EH Type) • **Setonic**

With Extended Handle

Cat. No.	Model	용량	Plunger Type	Needle		EA Price ₩
				od (Gauge)	길이	
S14-165-601	3031098	10 ul	SS	0.50mm (G25S)	50 mm	162,900
S14-165-605	3031099	25 ul				162,900
S14-165-609	3031100	100 ul				162,900
S14-165-613	3031108	250 ul				162,900
S14-165-617	3031103	25 ul	PTFE	0.50mm (G25S)	50 mm	162,900
S14-165-621	3031104	50 ul				162,900
S14-165-625	3031106	100 ul				162,900
S14-165-629	3031109	250 ul				162,900

Without Extended Handle

Cat. No.	Model	용량	Plunger Type	Needle		EA Price ₩
				od (Gauge)	길이	
S14-165-651	3031101	10 ul	SS	0.50mm (G25S)	50 mm	162,900
S14-165-655	3031102	25 ul				162,900
S14-165-659	3031105	50 ul				174,600
S14-165-663	3031107	100 ul				162,900
S14-165-667	3031110	250 ul				162,900


Gas Tight Syringe, TLL Luer Lock 가스 샘플 투입용 주사기

- 액상 및 기체상 시료에 모두 사용
- 테프론 팁이 유리 재질의 몸통과 빈틈없이 밀착하여 Tight한 Sealing 상태를 만들어줌
- 용도 : Sample 및 Aggressive Fluid의 Dosage
- Glass Barrel : Borosilicate 3.3 • Luer Lock Type
- 기타 제품은 www.labmart.co.kr 참조 • **Setonic**

Cat. No.	Model	용량	Plunger Type	Termination	EA Price ₩
S14-165-701	3010299	25 ul	PTFE	Luer Lock	84,000
S14-165-705	3010300	50 ul			84,000
S14-165-709	3010317	100 ul			84,000
S14-165-713	3010318	250 ul			88,700
S14-165-717	3010319	500 ul			86,400
S14-165-721	3010175	1 ml			71,700
S14-165-725	3010176	2.5 ml			79,800
S14-165-729	3010177	5.0 ml			79,200
S14-165-733	3010178	10.0 ml			88,900
S14-165-737	3010179	25.0 ml			202,700

Optional Metal Luer-Lock Needle

Cat. No.	Model	Gauge	길이	Type	PK Price ₩
S14-165-001	3050188	18G	50 mm	B 3ea	48,200
S14-165-003	3050266	19G	51 mm	A 3ea	34,900
S14-165-005	3050325			B 3ea	39,200
S14-165-007	3050383	D 3ea		65,100	
S14-165-009	3050446	A 1ea		20,000	
S14-165-011	3050441	20G	51 mm	B 3ea	39,200
S14-165-013	3050405			B 3ea	39,200
S14-165-015	3050410	21G		D 3ea	69,800
S14-165-017	3050104			B 3ea	39,200
S14-165-019	3050264	22G	51 mm	A 3ea	34,900
S14-165-021	3050379			AS 3ea	39,200
S14-165-023	3050380	D 3ea		65,100	
S14-165-025	3050262	23G		A 3ea	34,900
S14-165-027	3050323		AS 3ea	39,200	



Cat. No.	Model	Gauge	길이	Type	PK Price ₩
S14-165-029	3050355	23G	51 mm	B 3ea	39,200
S14-165-031	3050377			D 3ea	65,100
S14-165-033	3050198	B 3ea		23,300	
S14-165-035	3050260	A 3ea		34,900	
S14-165-037	3050357	26G	51 mm	B 3ea	39,200
S14-165-039	3050372			AS 3ea	39,200
S14-165-041	3050373			D 3ea	65,100
S14-165-043	3050223	27G		50 mm	A 3ea
S14-165-045	3050440		51 mm	B 3ea	39,200
S14-165-047	3050396	28G	25 mm	B 3ea	49,900
S14-165-049	3050222	29G	50 mm	A 3ea	99,800
S14-165-051	3050359	30G	51 mm	B 3ea	41,400
S14-165-053	3050360	32G		B 3ea	48,100
S14-165-055	3050361	33G		B 3ea	48,100