

Manual GC Syringe GC용 주사기

- 주용도 : Sample Dosage for Gas Chromatography, HPLC Valve, Autosampler
- Plunger : Stainless Steel
- Glass Barrel : Borosilicate 3.3
- Precision : $\langle \pm 1\%$ of the volume
- ILS

Fixed Needle Type ①

Cat. No.	Model	용량	포함된 Needle Tip			EA Price ₩
			Gauge	od	길이	
S15-204-202	2100301	5 ul	26S	0.47 mm	51 mm	92,900
S15-204-208	2100401	10 ul				221,100
S15-204-214	2100501	25 ul	22S	0.72 mm	51 mm	98,200
S15-204-220	2100601	50 ul				92,900
S15-204-226	2100701	100 ul	22			111,100
S15-204-232	2100801	250 ul				476,500
S15-204-238	2100901	500 ul				513,500

- Type A Needle 기본 포함

Removable Needle Type ②

Cat. No.	Model	용량	포함된 Needle Tip			EA Price ₩
			Gauge	od	길이	
S15-204-244	2100305	5 ul	26S	0.47 mm	51 mm	513,500
S15-204-250	2100405	10 ul				128,700
S15-204-256	2100505	25 ul	22S	0.72 mm	51 mm	442,700
S15-204-262	2100605	50 ul				548,600
S15-204-268	2100705	100 ul	22			150,500
S15-204-274	2100805	250 ul				638,200
S15-204-280	2100905	500 ul				638,200



Manual HPLC Syringe for Most Common Valve HPLC용 주사기, Common Valve용

- 주용도 : Sample Dosage for Gas Chromatography, HPLC Valve, Autosampler
- Plunger : Stainless Steel
- Glass Barrel : Borosilicate 3.3
- Precision : $\langle \pm 1\%$ of the volume
- ILS

Fixed Needle Type ③

Cat. No.	Model	용량	포함된 Needle Tip			EA Price ₩
			Gauge	od	길이	
S15-204-412	2101302	5 ul	22S	0.72 mm	51 mm	92,900
S15-204-418	2101402	10 ul				336,600
S15-204-424	2101502	25 ul				370,400
S15-204-430	2101602	50 ul				407,600
S15-204-436	2101702	100 ul				442,700
S15-204-442	2101802	250 ul	22			119,300
S15-204-448	2101902	500 ul				128,700

- Type B Needle 기본 포함

Removable Needle Type ④

Cat. No.	Model	용량	포함된 Needle Tip			EA Price ₩
			Gauge	od	길이	
S15-204-454	2101306	5 ul	22S	0.72 mm	51 mm	141,900
S15-204-460	2101406	10 ul				155,700
S15-204-466	2101506	25 ul				128,700
S15-204-472	2101606	50 ul				150,500
S15-204-478	2101706	100 ul				150,500
S15-204-484	2101806	250 ul	22			159,900
S15-204-500	2101906	500 ul				159,900

- Type B Needle 기본 포함

