

Multi-Purpose Labeling Tape 다목적 라벨 테이프

- 다목적으로 실험실에서 사용 가능한 라벨 테이프
- 거의 모든 실험기기에 부착 가능 : Glass, Plastic, Metal, Paper, Rubber 등
- Oilproof, Waterproof, Acid Resistant & Autoclared
- 페어낸 자리에 접착제가 남지 않으며, 접착후 깨끗하게 이동 가능
- Shamrock

길이 54.9 m, 3" Core, -20 ~ 110°C

Cat. No	Model	Size	색상	EA Price ₩
T01-218-201	ST-120-1	w 12.7 mm L 54.9 M	White	33,000
T01-218-203	ST-120-2		Yellow	33,000
T01-218-205	ST-120-3		Green	33,000
T01-218-207	ST-120-4		Red	33,000
T01-218-209	ST-120-5		Orange	33,000
T01-218-211	ST-120-6		Blue	33,000
T01-218-213	ST-120-7		Pink	33,000
T01-218-221	ST-340-1	w 19 mm L 54.9 M	White	45,600
T01-218-223	ST-340-2		Yellow	45,600
T01-218-225	ST-340-3		Green	45,600
T01-218-227	ST-340-4		Red	45,600
T01-218-229	ST-340-5		Orange	45,600
T01-218-231	ST-340-6		Blue	45,600
T01-218-233	ST-340-7		Pink	45,600
T01-218-241	ST-100-1	w 25.4 mm L 54.9 M	White	58,900
T01-218-243	ST-100-2		Yellow	58,900
T01-218-245	ST-100-3		Green	58,900
T01-218-247	ST-100-4		Red	58,900
T01-218-249	ST-100-5		Orange	58,900
T01-218-251	ST-100-6		Blue	58,900
T01-218-253	ST-100-7		Pink	58,900
T01-218-261	ST-110-1	w 38.1 mm L 54.9 M	White	83,000
T01-218-263	ST-110-2		Yellow	83,000
T01-218-265	ST-110-3		Green	83,000
T01-218-267	ST-110-4		Red	83,000
T01-218-269	ST-110-5		Orange	83,000
T01-218-271	ST-110-6		Blue	83,000
T01-218-273	ST-110-7		Pink	83,000

Dot Label Tape, 원형 라벨 테이프, -54 ~ 138°C

Cat. No	Model	Size	색상	EA Price ₩
T01-218-401	SDL-50R-1-REM	Ø12.7 mm, 1000/RL	White	25,200
T01-218-403	SDL-50R-2-REM		Yellow	25,200
T01-218-405	SDL-50R-3-REM		Green	25,200
T01-218-407	SDL-50R-4-REM		Red	25,200
T01-218-409	SDL-50R-5-REM		Orange	25,200
T01-218-411	SDL-50R-6-REM		Blue	25,200
T01-218-413	SDL-50R-7-REM		Pink	25,200


• Roll Type : 1000 Dot/Roll

Cryogenic Label Tape 초저온용 라벨 테이프, -196°C

- Cryopen-26으로 마감하기에 적합
- 사용 온도범위 : -196 ~ 90°C
- 색상 : White
- Shamrock

Roll Type, 500 Dot/Roll

Cat. No	Model	Size	EA Price ₩
T01-218-501	LCS-43R-D	Ø11.1 mm	53,100
T01-218-503	LCS-59R-D	12.7×22.9 mm	57,800
T01-218-507	LCS-759R-D	19.1×22.9 mm	70,100
T01-218-509	LCS-7512R-D	19.1×31.8 mm	76,800



길이 12.7 m, 1" Core, -20 ~ 110°C

Cat. No	Model	Size	색상	EA Price ₩
T01-218-301	ST-12-1	w 12.7 mm L 12.7 M	White	12,900
T01-218-303	ST-12-2		Yellow	12,900
T01-218-305	ST-12-3		Green	12,900
T01-218-307	ST-12-4		Red	12,900
T01-218-309	ST-12-5		Orange	12,900
T01-218-311	ST-12-6		Blue	12,900
T01-218-313	ST-12-7		Pink	12,900
T01-218-321	ST-34-1	w 19 mm L 12.7 M	White	17,100
T01-218-323	ST-34-2		Yellow	17,100
T01-218-325	ST-34-3		Green	17,100
T01-218-327	ST-34-4		Red	17,100
T01-218-329	ST-34-5		Orange	17,100
T01-218-331	ST-34-6		Blue	17,100
T01-218-333	ST-34-7		Pink	17,100
T01-218-341	ST-10-1	w 25.4 mm L 12.7 M	White	21,700
T01-218-343	ST-10-2		Yellow	21,700
T01-218-345	ST-10-3		Green	21,700
T01-218-347	ST-10-4		Red	21,700
T01-218-349	ST-10-5		Orange	21,700
T01-218-351	ST-10-6		Blue	21,700
T01-218-353	ST-10-7		Pink	21,700
T01-218-361	ST-11-1	w 38.1 mm L 12.7 M	White	30,700
T01-218-363	ST-11-2		Yellow	30,700
T01-218-365	ST-11-3		Green	30,700
T01-218-367	ST-11-4		Red	30,700
T01-218-369	ST-11-5		Orange	30,700
T01-218-371	ST-11-6		Blue	30,700
T01-218-373	ST-11-7		Pink	30,700

Cat. No	Model	Size	색상	EA Price ₩
T01-218-421	SDL-100R-1-REM	Ø25.4 mm, 1000/RL	White	33,600
T01-218-423	SDL-100R-2-REM		Yellow	33,600
T01-218-425	SDL-100R-3-REM		Green	33,600
T01-218-427	SDL-100R-4-REM		Red	33,600
T01-218-429	SDL-100R-5-REM		Orange	33,600
T01-218-431	SDL-100R-6-REM		Blue	33,600
T01-218-433	SDL-100R-7-REM		Pink	33,600

• Roll Type : 1000 Dot/Roll



Sheet Type, 500 ea/Sheet

Cat. No	Model	Size	EA Price ₩
T01-218-521	LCS-43D	Ø11.1 mm	53,100
T01-218-523	LCS-59D	12.7×22.9 mm	57,800
T01-218-527	LCS-759D	19.1×22.9 mm	70,100
T01-218-529	LCS-7512D	19.1×31.8 mm	76,800