

Adapter, Glass Thermometer Holder 유리 온도계 어댑터

Cat. No.	Model	Thread	Cone	EA Price ₩
A01-83-608	1311414S	GL 14	14/23	18,700
A01-83-612	1311419S		19/26	21,800
A01-83-616	1311424S		24/29	26,500
A01-83-620	1311429S		29/32	26,500

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- GL14 Open top Screw Cap & PTFE Faced Silicon Rubber Septa 포함
- Bore : 5.5 ~ 6.5 mm
- Lenz

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
B08-83-852	1321014	Open-top Screw Cap, PBT, GL 14	3,400
J02-83-002	1331406	PTFE Faced Silicon Rubber Septa	5,300



Adapters, Thermometer Inlet 유리 온도계 어댑터

Cat. No.	Model	Joint	높이	Thread	EA Price ₩
A01-132-810	CG-1042-14	14/20	50mm	#7 (od4-7mm용)	60,200
A01-132-813	CG-1042-19	19/22	65mm		62,700
A01-132-816	CG-1042-01	24/40	80mm		61,400
A01-132-819	CG-1042-02	29/42	85mm		68,900
R02-132-956	CG-1042-45	45/50	95mm	#15 (od12-15mm용)	117,500
R02-132-984	CG-1042-46	45/50	95mm		121,800

- Compression Cap에 Tube 연결을 위한 Hole이 있음
- od 4-7mm Tube(#7 Cap) 또는 od 12-15mm Tube(#15 Cap)와 같이 사용하기에 맞맞음
- Chem-Thread Compression Cap 및 Viton O-ring 포함
- Chemglass

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
A01-132-843	CG-352-01	#7 Cap for 4~7mm od Tube	7,700
A01-132-852	CG-352-03	#15 Cap for 12~15mm od Tube	16,000
A01-132-870	CG-305-107	Viton O-Ring, #107, 12/pk	19,300
A01-132-873	CG-309-107	Perfluoro O-Ring, #107, Black	69,000
A01-132-876	CG-305-206	Viton O-Ring, #206, 6/pk	22,300
A01-132-879	CG-309-206	Perfluoro O-Ring, #206, Black	87,700



Bushing Type Thermometer Adapter 부싱 타입 어댑터, LukeGL®

Cat. No.	Model	적용 외경	Cone	EA Price ₩
A01-132-802	360-01-24	∅6 ~ 7 mm	24/40	32,000
A01-132-804	360-03-24	∅12.5 ~ 14 mm		50,600

- Borosilicate Glass 3.3
- 부싱타입 온도계 어댑터 (360-01-24)
- 상부 Fitting 재질 : Nylon (360-01-24), PTFE (360-03-24)
- 구성 : Glass Body, 상부 Fitting, Viton O-ring
- LukeGL®

