

## I Laboratory Blender 실험실용 분쇄기, 1 L

• Waring

### Model : A type - 속도 고정형

Cat. No.	Model	Container	Motor & Blade	EA Price ₩
B06-74-006	700S	SUS, 1 L	0.40 Hp Standard	980,700
B06-74-012	700G	Glass, 1.2 L		806,000

- Two-piece Vinyl and Styrene Lid 포함
- Epoxy-Coated Motor Housing
- 속도 고정형 : 22,000 rpm
- Power : 120V, 60 Hz

### Model : B type - 속도 가변형

Cat. No.	Model	Container	Motor & Blade	EA Price ₩
B06-74-018	7011S	SUS	0.40 Hp Standard	1,149,400
B06-74-024	7011G	Glass		1,047,000
B06-74-030	7011HS	SUS	0.70 Hp Heavy Duty	1,360,200
B06-74-036	7011HG	Glass		1,257,800

- Two-piece Vinyl and Styrene Lid 포함
- Epoxy-Coated Motor Housing
- 속도 가변형 : Two Speed (High 22,000 rpm, Low 18,000 rpm)
- Power : 120V, 60 Hz

### Model : C type - 속도 가변형

Cat. No.	Model	Container	Motor & Blade	EA Price ₩
B06-74-042	LB10S	0.40 Hp SUS	0.40 Hp Standard	2,097,600
B06-74-048	LB10G	0.40 Hp Glass		1,975,900

- Two-piece Vinyl and Styrene Lid 포함
- 속도 가변형 : Variable Speed (1,000 ~ 22,000 rpm, with 외부 조절 장치)
- Power : 120V, 60 Hz

### Model : D type - 속도 가변형

Cat. No.	Model	Container	Motor & Blade	EA Price ₩
B06-74-066	7010S	SUS	0.40 Hp Standard	1,360,200
B06-74-072	7010G	Glass		1,219,300
B06-74-078	7010HS	SUS	0.70 Hp Heavy Duty	1,622,900
B06-74-084	7010HG	Glass		1,520,500

- Two-piece Vinyl and Styrene Lid 포함 (Epoxy-Coated Motor Housing)
- 속도 가변형 : Two Speed with 3 min Timer (High 22,000, Low 18,000)
- Power : 120V, 60 Hz

A type - 속도 고정형, SUS




A type - 속도 고정형, Glass



B type - 속도 가변형, SUS



B type - 속도 가변형, Glass



C type - 속도 가변형, SUS



C type - 속도 가변형, Glass



D type - 속도 가변형, SUS



D type - 속도 가변형, Glass

