

K-type Thermometer Probe / Thermocouple K 타입 온도 센서 / 열전대

Cat. No	Model	Description	EA Price ₩
T22-147-792	4014	Triple Purpose Probe, -50 ~ 700°C	118,300
T22-147-799	4028	Instand Probe (Beaded Probe), Teflon Coating	54,700


- K-Type 온도 센서 : Nicar / NiAl
- Triple Purpose : Liquid, Air /Gas, Semi-Solid용
- Cole-Parmer / Traceable®

센서 사양

Cat. No	Model	온도 범위	Stem Size	전체 길이	Cable
T22-147-792	4014	-50 ~ 700°C	∅3.3×L159 mm	248 mm	1 M (길이)
T22-147-799	4028	-40 ~ 250°C	-	1.2 M	∅ 1.5 mm



4014



4028

K-type Thermometer Probe K 타입 열전대 ①

Cat. No	Model	Length	Stem	EA Price ₩
T22-132-754	L349820	200 mm	∅6.4 mm	61,500
T22-132-760	L349830	300 mm		61,500
T22-132-766	L349840	400 mm		65,600
T22-132-772	L349850	500 mm		65,600
T22-132-778	L349860	600 mm		73,200

- 재질 : Stainless Steel
- K-Type Connector를 가지고 있는 온도센서
- LK Labkorea

Thermocouple Extension Cord ②

Cat. No	Model	Length	EA Price ₩
T22-132-882	L349902	2 M Coil	36,200

Thermocouple Clamp ③

Cat. No	Model	Description	EA Price ₩
T22-132-955	3499-H-10	1.6 or 3.2 or 6.4 mm Rod 용	145,400


- Chemglass



①



②



③

K-type Thermocouple Probe K 타입 온도 센서

Cat. No	Model	Size	Ferrule	Cable	EA Price ₩
T22-450-101	600.106.1	∅6×L100 mm	35 mm	1 M	100,500
T22-450-105	600.206.1	∅6×L200 mm	40 mm		109,300
T22-450-109	600.101.1	∅7×L100 mm	35 mm		107,000
T22-450-113	600.201.1	∅7×L200 mm	40 mm		116,000
T22-450-117	600.307.1	∅7×L300 mm			118,200
T22-450-121	600.208.1	∅8×L200 mm			111,700
T22-450-125	600.301.1	∅8×L300 mm			125,800
T22-450-129	600.401.1	∅8×L400 mm	50 mm		134,100
T22-450-133	600.501.1	∅9×L500 mm	60 mm		153,300
T22-450-137	600.601.1	∅9×L600 mm			168,300
T22-450-141	600.701.1	∅9×L700 mm			194,900
T22-450-145	600.901.1	∅9×L900 mm	45 mm		198,800

- 재질 : PTFE, -180 ~ 280°C 내열
- Supplied with 1m Cable
- Cowie

