

Dounce Tissue Grinder, All Glass, Kimble®

다운스 티슈 그라인더, 2×Pestle

Cat. No	Model	용량	Pestle Clearance		EA Price ₩
			Large	Small	
T24-137-011	885300-0000	0.5 ml	0.0025-0.0055	0.0005-0.0025	320,100
T24-137-015	885300-0001	1 ml			320,100
T24-137-019	885300-0002	2 ml	0.0030-0.0050		367,100
T24-137-023	885300-0007	7 ml	0.0028-0.0047	0.0008-0.0022	324,700
T24-137-027	885300-0015	15 ml	0.0035-0.0065	0.0010-0.0030	364,500
T24-137-031	885300-0040	40 ml	0.0030-0.0060		421,300
T24-137-035	885300-0100	100 ml	0.0020-0.0100	0.0005-0.0055	906,400

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Pestle Clearance : Tube 와 Pestle 사이의 공간
- DWK Lifescience / Kimble®



DUALL® Tissue Grinder, Conical Shape, Kimble®

듀얼 코니칼 티슈 그라인더

- DWK Lifescience / Kimble®

Micro DUALL® All Glass ①

Cat. No	Model	Working Capacity	Pestle Shaft	PK Price ₩
T24-137-039	885470-0000	0.2 ml	od3 mm	6ea 755,900

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 제품 전체가 유리로 제작됨
- 버퍼액 속에 있는 적은 양의 Tissue 분쇄에 사용
- PTFE Bearing은 Aerosol이 생기는 것과 위치가 어긋나는 것을 방지



Micro DUALL® All Glass with Screwcap ②

Cat. No	Model	Working Capacity	Pestle Shaft	EA Price ₩
T24-137-139	885490-0017	0.3 ml	od6 mm	261,900
T24-137-143	885490-0019	1 ml		277,000
T24-137-147	885490-0020	3 ml	od9 mm	320,200
T24-137-151	885490-0021	5 ml		334,700

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 원심 분리, 결합, 냉동 건조 등에도 사용 가능
- Heavy Wall 타입으로 손상이 적음 • Pestle과 Tube 사이의 공간 : 0.010 ~ 0.015 mm



DUALL® with All Glass ③

Cat. No	Model	Working Capacity	Pestle Shaft	EA Price ₩
T24-137-043	885450-0020	1 ml	od6 mm	313,700
T24-137-047	885450-0021	3 ml		313,700
T24-137-051	885450-0022	5 ml	od8 mm	316,300
T24-137-055	885450-0023	15 ml	od10 mm	361,000
T24-137-059	885450-0024	30 ml		384,600
T24-137-063	885450-0025	50 ml	od16 mm	488,200

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Grinding 효율을 최상으로 극대화
- 근육, 심장, 폐 등 결합 조직과 사용하기에 적합
- Pestle과 Tube 사이의 공간 : 0.010 ~ 0.015 mm

DUALL® with Glass Pestle & Kimcote® Glass Tube ④

Cat. No	Model	Working Capacity	Pestle Shaft	EA Price ₩
T24-137-067	885460-0020	1 ml	od6 mm	329,600
T24-137-071	885460-0021	3 ml		329,600
T24-137-075	885460-0022	5 ml	od8 mm	329,600
T24-137-079	885460-0023	15 ml	od10 mm	402,700
T24-137-083	885460-0024	30 ml		497,600
T24-137-087	885460-0025	50 ml	od16 mm	612,900

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Glass Tube는 플라스틱으로 코팅되어 있어 Tube가 깨지거나 금이 가도 균질화 작업 지속 가능
- Pestle과 Tube 사이의 공간 : 0.010 ~ 0.015 mm

