

Mercury Tong 수은 집게

Cat. No	Model	Grip Capacity	길이	EA Price ₩
T19-139-331	1488-22	25 ~ 80 mm	220 mm	24,900
T19-139-333	1488-25	25 ~ 90 mm	250 mm	27,300

- 재질 : Stainless Steel 304
- LK Labkorea



Laboratory Tong 실험실용 집게

Cat. No	Model	Grip Capacity	길이	EA Price ₩
T19-139-351	1494-40	10 ~ 130 mm	400 mm	24,900


- 재질 : Stainless Steel 304
- LK Labkorea



Utility Straight Tong 다용도 직선형 집게

Cat. No	Model	Opening Size	길이	EA Price ₩
T19-139-141	1454-20	20 mm	150 mm	12,400
T19-139-143	1454-22	20 mm	200 mm	14,200
T19-139-145	1454-25	25 mm	250 mm	16,100
T19-139-147	1454-30	25 mm	300 mm	17,500
T19-139-149	1454-40	40 mm	400 mm	22,400
T19-139-151	1454-50	40 mm	500 mm	28,700
T19-139-153	1454-60	40 mm	600 mm	36,000

- 재질 : Stainless Steel 304
- LK Labkorea



Sterilizing Tong 멸균기용 집게

Cat. No	Model	Grip Capacity	길이	EA Price ₩
T19-139-103	1490-28	100 ~ 200 mm	280 mm	17,500

- 재질 : Stainless Steel 304
- LK Labkorea



Tube Tong 튜브 집게

Cat. No	Model	Grip Capacity	EA Price ₩
T19-284-122	4160	Width of Mouth 59 mm	52,300

- 재질 : 니켈 도금
- 아연 다이캐스팅
- Jaw 부분은 비닐로 코팅되어 있어 미끄러짐 방지하며, 열과 화학물질에 강함
- Usbeck



Test Tube Holder with Finger 시험관 집게

Cat. No	Model	Grip Capacity	길이	EA Price ₩
T19-139-391	1548-13	8 ~ 90 mm	130 mm	4,200
T19-139-393	1548-15	8 ~ 70 mm	150 mm	4,900
T19-139-395	1548-18	8 ~ 60 mm	180 mm	6,300

- 재질 : Stainless Steel 304
- LK Labkorea

