



**Glass Rectangular Tray, Simax®**  
**유리 직사각 트레이**

Cat. No	Model	용량	Size (mm)		EA Price ₩
			w×d	h	
T09-985-311	31121202	2 L	200×100	100	53,900
T09-985-313	31121203	1 L	200×100	50	41,100
T09-985-315	31121204	2.5 L	250×165	60	53,000
T09-985-317	31121205	5 L	250×165	120	80,100

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 사용온도 : -40 ~ 300°C
- KAVALIERGLASS



**Glass Round Tray with Lid, Simax®**  
**유리 원형 트레이**

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 사용온도 : -40°C ~ 300°C
- Autoclave 가능하며, Microwave Oven에서도 사용 가능함
- KAVALIERGLASS



**Low Form, 낮은 유리 원형 트레이**

Cat. No	Model	용량	Size	EA Price ₩
T09-985-101	6056 / 6066	700 ml	Ø160×h70 mm	13,900
T09-985-103	6256 / 6266	1 L	Ø182×h81 mm	15,700
T09-985-105	6106 / 6116	1.5 L	Ø205×h93 mm	17,200
T09-985-107	6026 / 6036	2 L	Ø228×h102 mm	20,600
T09-985-109	6176 / 6186	2.5 L	Ø249×h113 mm	24,200
T09-985-111	6226 / 6236	3.5 L	Ø268×h118 mm	28,200

**Deep Form, 깊은 유리 원형 트레이**

Cat. No	Model	용량	Size	EA Price ₩
T09-985-121	6666 / 6676	1.5 L	Ø188×h100 mm	18,700
T09-985-123	6686 / 6696	1.75 L	Ø208×h113 mm	21,400
T09-985-125	6706 / 6716	2.8 L	Ø228×h125 mm	26,300

**Glass Oval Tray, Simax®**  
**유리 타원형 트레이**

Cat. No	Model	용량	Size (w×d×h)	EA Price ₩
T09-985-131	7127	2.2 L	300×210×62 mm	15,700
T09-985-133	7126	3 L	350×240×62 mm	19,800
T09-985-135	7276	4 L	400×270×62 mm	25,800

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 사용온도 : -40°C ~ 300°C
- Autoclave 가능하며, Microwave Oven에서도 사용 가능함
- KAVALIERGLASS

