

Plastic Utility Tray / Sterilizing Pan 플라스틱 트레이

| Cat. No | Model | w (mm) | d (mm) | h (mm) | EA Price ₩ |
|------------|---------|--------|--------|--------|------------|
| T09-46-931 | T242000 | 320 | 260 | 70 | 10,800 |
| T09-46-933 | T242010 | | | 100 | 13,600 |
| T09-46-935 | T240070 | 360 | 310 | 130 | 19,500 |
| T09-46-937 | T240090 | 540 | 435 | | 26,700 |

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Sterilizing Pan으로 사용 할 수 있음
- 손잡이가 달려 있어 이동이 편리함
- Tarsons



PP Laboratory Tray 플라스틱 트레이

| Cat. No | Model | w (mm) | d (mm) | h (mm) | EA Price ₩ |
|------------|---------|--------|--------|--------|------------|
| T09-53-258 | 1050-94 | 225 | 180 | 45 | 19,400 |
| T09-53-264 | 1051-94 | 300 | 240 | 70 | 32,900 |
| T09-53-270 | 1052-94 | 370 | 310 | 75 | 47,100 |
| T09-53-276 | 1053-94 | 510 | 410 | 120 | 69,600 |
| T09-53-282 | 1054-94 | 625 | 530 | 140 | 120,700 |
| T09-53-288 | 1055-94 | 840 | 645 | 160 | 175,900 |

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Brand GmBH



PVC Tray PVC 다용도 트레이

| Cat. No | Model | Size (w×h, mm) | | Height | EA Price ₩ |
|-------------|----------|----------------|---------|--------|------------|
| | | Inner | Outer | | |
| T09-100-011 | 42021013 | 140×190 | 180×230 | 40 mm | 21,900 |
| T09-100-015 | 42021018 | 190×260 | 240×320 | 50 mm | 33,100 |
| T09-100-019 | 42021024 | 260×320 | 320×370 | 50 mm | 51,600 |
| T09-100-023 | 42021030 | 310×410 | 370×470 | 60 mm | 85,400 |
| T09-100-027 | 42021040 | 420×510 | 480×580 | 60 mm | 162,700 |
| T09-100-031 | 42021050 | 510×610 | 590×680 | 70 mm | 231,000 |

- 재질 : PVC, -30 ~ 70°C 내열
- 사진 현상용이나 일반 Tray의 용도로 사용됨
- Acid나 Alkalis에 강한 내성을 가짐
- Wenk Labtec



Compartment Tray 사각 칸막이 트레이

| Cat. No | Model | 칸수 | 칸 규격, w×d | EA Price ₩ |
|------------|--------|----|-----------------|------------|
| T09-20-659 | 768510 | 12 | 90×90 mm (12칸) | 101,900 |
| T09-20-663 | 768505 | 5 | 385×95 mm (1칸) | 101,900 |
| | | | 190×120 mm (2칸) | |
| | | | 190×90 mm (2칸) | |

- 재질 : PVC, -30 ~ 70°C 내열
- 전체 Size : w400×d300×h65 mm
- Brand GmBH



T09-20-663