

Laboratory Blender 실험실용 분쇄기

Cat. No.	Model	용량	Speed	EA Price ₩
B06-74-310	WB-1	150 ml	25000 rpm	1,156,600
B06-74-410	WDL-1	200 ml		3,177,300

- 곡물, 약재, 광물, 목재 등 연질의 시료 분쇄용
- Container 재질 : Stainless Steel
- Wonder Blender

Cat. No.	B06-74-310	B06-74-410
Model	WB-1	WDL-1
용량	150 ml	200 ml
Rotatin Speed	25000 rpm	
Power Consumption	700 W	1 kW
Size	Ø134×h270 mm	Ø162×h298 mm
Weight	3.6 Kg	4.6 Kg
Accessory	Minute Crush Lid. Brush	Brush
Power	220 V, 50/60 Hz	110 V, 50/60 Hz



Laboratory Crusher / Blender 실험실용 분쇄기

Cat. No.	Model	용량	Max. Speed	Minute Lid	EA Price ₩
B06-74-510	WC-3 ①	150 ml	28000 rpm	X	2,156,000
B06-74-520	WC-3L ②			O	2,451,000

- 28000 rpm의 Motor를 장착한 High Performance Crusher
- Rock, Dry Crude Drug, Medicinal Herbs, Tablet, Mineral 및 Wood Chip의 분쇄용
- Stainless Steel Container와 Motor Part (Main Body)가 분리되어짐
- 냉각 기능이 내장된 용기를 사용 할 수 있어 열에 민감한 재료의 파쇄에 적합함 (Option : PN-L15C)
- Wonder Blender

Cat. No.	B06-74-510	B06-74-520
Model	WC-3	WC-3L
Type	with Standard Lid	with Minute Crush Lid
용량	150 ml	
Rotatin Speed	0 ~ 28000 rpm (No Load) / Variable	
Power Consumption	1.2 kW	
작동 시간	1 min	
Size	용기	Ø108×h51 mm
	전체	Ø170×h330 mm
Weight	5 Kg	
Power	AC100 V, 50/60 Hz	

Container Set

Cat. No.	Model	용량	Description	EA Price ₩
B06-74-530	PN-L15 ③	150 ml	SUS Container with Standard Lid	1,248,200
B06-74-531	PN-L15C	150 ml	SUS Cooling Container with Standard Lid	1,621,100
B06-74-535	PN-L16 ④	75 ml	SUS Container with Standard Lid	1,134,800
B06-74-537	PN-W03		Crushed Lid	527,100

