

## PFA Tray PFA 테프론 트레이

Cat. No	Model	외부 Size (mm)	내부 Size (mm)	EA Price ₩
T09-64-122	F1194-01	w242×d171×h40	w205×d134×h35	343,700
T09-64-126	F1194-02	w304×d216×h43	w264×d175×h39	485,500
T09-64-130	F1194-03	w350×d252×h46	w309×d210×h41	823,700

- 재질 : PFA
- Flon Industry



## PFA Small Tray PFA 테프론 소형 트레이

Cat. No	Model	외부 Size (mm)	내부 Size (mm)	EA Price ₩
T09-64-222	F1193-01	w106×d78×h27	w98×d70×h24	83,900
T09-64-226	F1193-02	w119×d88×h27	w112×d80×h24	94,000
T09-64-230	F1193-03	w137×d102×h27	w128×d94×h23	132,900

- 재질 : PFA
- Flon Industry



## Conductive Tray 정전기 방지용 트레이

Cat. No	Model	w (mm)	d (mm)	h (mm)	EA Price ₩
T09-143-129	CA421C	485	355	45	16,100
T09-143-136	CA422C	355	255	45	14,500
T09-143-143	CA423C	280	215	40	13,000
T09-143-150	CA424C	200	150	40	9,600

- 재질 : Conductive PS, Black
- 전기저항 : 104~106 Ω/Sq
- LK Labkorea



## PP Large Tray PP 대형 트레이와 커버

Cat. No	Model	w (mm)	d (mm)	h (mm)	EA Price ₩
T09-517-880	39512	380	500	125	16,100
T09-517-888	39518	390	520	180	21,500

- 재질 : Polypropylene (PP)
- 포개어서 보관 할 수 있음
- 커버는 모든 사이즈와 호환 가능하며 별도 구매해야함
- Rayhope®

### Optional Lid

Cat. No	Model	Description	EA Price ₩
T09-517-950	39520	w392×d522×h20 mm	11,700

