



## Centrifuge Tube, ASTM ASTM 원심관

Cat. No.	Model	용량	ASTM	Max. Length	EA Price ₩
T16-52-905	362038 ①	100 ml	D91	203 mm	78,300
T16-52-903	362338 ②		D96	167 mm	81,000
T16-52-900	362138 ③			160 mm	127,000

- Borosilicate Glass 3.3
- Max. RCF : 700 (at 20 °C for 20 min)
- ASTM 규격에 따름
- Brand GmBH

### 용량별 눈금

Cat. No.	T16-52-905		T16-52-903		T16-52-900	
	362038		362338		362138	
단위	용량, ml	눈금, ml	용량, ml	눈금, ml	용량, ml	눈금, ml
	0 ~ 0.5	0.05	0 ~ 0.5	0.05	0 ~ 1.5	0.10
	0.5 ~ 2	0.10	0.5 ~ 2	0.10	1.5 ~ 3	0.50
	2 ~ 3	0.20	2 ~ 3	0.20	3 ~ 5	0.50
	3 ~ 5	0.50	3 ~ 5	0.50	5 ~ 10	1
	5 ~ 10	1	5 ~ 10	1	10 ~ 25	5
	10 ~ 25	5	10 ~ 25	5	25 ~ 100	25
	25 ~ 100	25	25 ~ 100	25	-	-



## Conical Bottom Centrifuge Tube 15 ml 코니칼 원심관

Cat. No.	Model	용량	PK Price ₩
T16-52-918	778013	15 ml	10ea 23,700

- Borosilicate Glass 3.3, USP Type I & ASTM E438 Type, Class B
- Size: Ø17×h113 mm, 두께 0.1 mm
- Max. RCF : 3500 (at 20 °C for 20 min)
- Brand GmBH



## Conical Bottom Centrifuge Tube 유리 코니칼 원심관, LukeGL®

Cat. No.	Model	용량	Size (od×h)	PK Price ₩
T16-256-202	CT-CB12T	12 ml	16.4×120 mm	5ea 28,000
T16-256-209	CT-CB15T	15 ml	16.4×135 mm	5ea 28,000

- Borosilicate Glass 3.3
- 멸균가능한 Teflon/Rubber Lined Screw Cap 일체형
- Screw Cap Size : 15-415
- 0.1 ml 눈금
- LukeGL®



## Round Bottom Centrifuge Tube 라운드 바텀 원심관, LukeGL®

Cat. No.	Model	용량	Size (od×h)	PK Price ₩
T16-256-655	CT-RB25T	25 ml	24×110 mm	5ea 27,100
T16-256-662	CT-RB40T	40 ml	28×110 mm	5ea 36,900

- Borosilicate Glass 3.3
- 멸균 가능한 Teflon/Rubber Lined Screw Cap 일체형
- Screw Cap Size : 24-400
- LukeGL®

