

Graduated Centrifuge Tube, Kimble® 눈금형 원심관


- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Calibrated to Contain (TC)
- Max. RCF : 2,980 xg
- DWK Lifescience / Kimble®

기본형 ①

Cat. No.	Model	용량	Size	PK Price ₩
T16-137-151	45165-10	10 ml	Ø17×L108 mm	12ea 461,900
T16-137-155	45165-15	15 ml	Ø17×L118 mm	12ea 318,600
T16-137-159	45165-50	50 ml	Ø29×L133 mm	12ea 533,300

Heavy-duty Type ②

Cat. No.	Model	용량	Size	PK Price ₩
T16-137-167	45199-12	12 ml	Ø17×L118 mm	12ea 422,700



Graduated Centrifuge Tube with Screw Cap, Kimble® 스크류 캡 원심관 ③

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
T16-137-179	45166-15	15 ml	17×130 mm	15-415	12ea 498,400
T16-137-183	45166-50	50 ml	29×147 mm	24-410	12ea 760,700
T16-137-187	45246-100	100 ml	37×203 mm	28-410	12ea 1,112,400

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Calibrated to Contain (TC)
- Black Phenolic white Rubber Liner Cap
- DWK Lifescience / Kimble®



High Strength Centrifuge Tube, Kimble® 고강도 원심관

- 재질 : Borosilicate Glass (only Tube) • ASTM E438, Type I, Class B
- Max. Temp. : up to 300°C (Only Tube)
- Max. RCF : 13,100 xg (with Rubber Sleeve in a 50 ml Rotor Cavity)
- 일반 유리보다 강한 강도와 내화학성이 있음
- DWK Lifescience / Kimble®

Plain Type ①

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩
T16-137-207	45500-15	15 ml	18×102 mm	6ea 136,700
T16-137-211	45500-30	30 ml	24×106 mm	6ea 219,200


Screw Thread Type ②

Cat. No.	Model	용량	Size (Ø×h)	Cap	PK Price ₩
T16-137-199	45600-15	15 ml	18×102 mm	20-400	6ea 290,400
T16-137-203	45600-30	30 ml	24×106 mm	24-400	6ea 294,300

- Black Phenolic PTFE-faced Rubber Liner Cap 포함

Rubber Adapter Sleeve ③

Cat. No.	Model	Description	PK Price ₩
T16-137-722	45550-15	15 ml Centrifuge Tube 용	2ea 67,300
T16-137-725	45550-30	30 ml Centrifuge Tube 용	2ea 76,500



Heavy-duty Round Bottom Centrifuge Tube with Cap 스크류 캡 헤비윌 원심관, Kimble®

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩
T16-137-745	45212-35	35 ml	29×100 mm	12ea 301,200
T16-137-748	45212-50	50 ml	29×123 mm	12ea 425,500

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Black Phenolic Rubber Liner Cap 포함
- DWK Lifescience / Kimble®

