

Alumina Boat 알루미나 보트

- 사용온도 : 1,600°C (Max. Temp. : 1,750°C)
- 내약품성이 뛰어나 화학적 침식에 높은 저항을 갖고 있음
- 산과 알칼리, 유기용매에 거의 영향을 받지 않음
- 일반 금속보다 내마모성이 15 ~ 20배 높음
- 반도체, LCD, 일반 산업용
- LK Labkorea

General Type


Cat. No.	Model	Size	EA Price ₩
B07-110-010	LCB003	w10×d50×h9 mm	7,100
B07-110-020	LCB005	w13×d70×h10 mm	7,800
B07-110-040	LCB015	w20×d100×h13 mm	9,200
B07-110-070	LCB050	w40×d80×h32 mm	18,200
B07-110-075	LCB055	w33×d113×h20 mm	31,400
B07-110-080	LCB060	w30×d136×h21 mm	31,400
B07-110-090	LCB075	w44×d160×h31 mm	48,200

U-Shape Type

Cat. No.	Model	Size	EA Price ₩
B07-110-110	ULCB003	w24×d32×h17 mm	7,100
B07-110-120	ULCB005	w45×d90×h22 mm	35,200
B07-110-140	ULCB015	w60×d180×h30 mm	52,700

Cover for U-Shape Boat

Cat. No.	Model	Description	EA Price ₩
B07-110-210	ULCB003C	Cover for ULCB003	9,100
B07-110-220	ULCB005C	Cover for ULCB005	30,200
B07-110-230	ULCB015C	Cover for ULCB015	48,200



Porcelain Boat Combustion 자제 연소용 보트

Cat. No.	Model	Size	EA Price ₩	PK Price ₩
B07-128-311	214/2	w75×d13×h8 mm	1,300	50ea 60,100
B07-128-321	214/4	w85×d13×h8 mm	1,300	50ea 60,100
B07-128-331	214/8	w90×d13×h8 mm	1,300	50ea 60,100
B07-128-341	214/3	w105×d14×h9 mm	1,300	50ea 60,100

- 재질 : Porcelain, 1,000°C 내열
- Jipo / LabCermic®



Quartz Boat 석영 보트

Cat. No.	Model	용량	Size	EA Price ₩
B07-60-009	LQB102	2 ml	w15×L48×h8 mm	33,900
B07-60-018	LQB103	3 ml	w16×L77×h10 mm	45,600
B07-60-027	LQB104	5 ml	w20×L77×h11 mm	54,700
B07-60-036	LQB108	8 ml	w20×L92×h11 mm	65,000
B07-60-045	LQB110	10 ml	w20×L102×h11 mm	78,000
B07-60-053	LQB112	15 ml	w27×L105×h11 mm	108,400

- 사용온도 : 1,100°C (Max.Temp. : 1,250°C)
- 비철금속 및 고온 물질 용해시 사용
- LK Labkorea

