

## Heavy Duty Glass Centrifuge Bottle, Kimble® 헤비월 유리 원심분리기용 바틀

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
B08-137-152	14720-200	200 ml	60×145 mm	105,300

- KIMAX® Borosilicate Glass 3.3, USP Type I, ASTM E438 Type I, Class A
- Small Bottom Area : Better Concentration, Max.RCF : 3550
- Pulp / Vinyl Coated Liner를 가진 Black Closure, GPI : 38-400
- DWK Lifescience / Kimble®



## Oak Ridge Centrifuge Tube, PP 오크리지 고속 원심관

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩
T16-46-052	T541020	10 ml	16×80.3 mm	12ea 20,600
T16-46-060	T541030	30 ml	25.5×96.3 mm	12ea 28,300
T16-46-068	T541040	50 ml	28×110 mm	12ea 38,200
T16-46-069	T541050	70 ml	37.5×102.2 mm	12ea 55,800

- 재질 : 반투명 (PP Body + PP Cap), -10 ~ 135°C 내열, Autoclavable
- Max.RCF : 50,000 xg
- Tarsons



## Oak Ridge Centrifuge Tube, FEP 테프론 오크리지 고속 원심관

Cat. No.	Model	Size Ø×h	Max.RPM	EA Price ₩
T16-64-058	F1122-001	16×80 mm	21,000	79,800
T16-64-061	F1122-002	25.5×92 mm	20,000	112,900
T16-64-064	F1122-003	28.5×104 mm	19,000	135,800

- 재질 : FEP 테프론, FEP Screw Cap
- -100°C ~ 150°C
- Flon Industry



## PFA Conical Centrifuge Tube PFA 테프론 코니칼 원심관

Cat. No.	Model	용량	Size	EA Price ₩
T16-64-801	F1182-001	15 ml	Ø16×h125 mm	252,800
T16-64-803	F1182-002	50 ml	Ø29×h115 mm	341,900
T16-64-805	F1182-003	55 ml	Ø29×h119 mm	299,200

- 재질 : PFA, -196 ~ 250°C 내열, Autoclavable
- Max. RCF : 50,000G
- Flon Industry



## Self-standing Conical Tube, 50ml 50ml 셀프 스텐딩 코니칼 튜브

Cat. No.	Model	Capacity	Size (Ø×h)	내포장	CS Price ₩
T16-35-051	CT-S050	50 ml	24×110 mm	25ea/pk 20pk	177,100

- 재질 : PP Tube (우수한 내구력) / PE Cap
- Sterilization : DNase, RNase Free, Non-Pyrogenic
- 지퍼백 내포장으로 재밀폐 가능 (오염 방지)
- RCF Rating : Max. 16000 xg
- LK Labkorea

