

PP Centrifuge Tube, Round Bottom 플라스틱 원심관, 뚜껑 별매

Cat. No.	Model	용량	Size (Ø×h)	PK Price ₩
T16-59-103	115342	10 ml	Ø16×100 mm	250ea 198,200
T16-59-106	115346	26 ml	Ø24×90 mm	50ea 113,000
T16-59-109	115348	30 ml	Ø20×100 mm	50ea 134,900
T16-59-112	115350	48 ml	Ø30×100 mm	25ea 69,600
T16-59-115	115352	75 ml	Ø35×100 mm	20ea 58,200
T16-59-118	115354	100 ml	Ø40×120 mm	20ea 66,300
T16-59-121	115356	160 ml	Ø45×120 mm	10ea 47,600

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Max. RCF : 4500 xg • 입구에 Rim이 있어 Stopper와 같이 사용하기에 적합
- Brand GmbH

PE Stopper for Centrifuge Tube

Cat. No.	Model	규격	PK Price ₩
T16-59-133	115360	10 ml (115342)용	100ea 10,100
T16-59-136	115366	30 ml (115348)용	50ea 10,500
T16-59-139	115368	26 ml (115346)용	50ea 8,600
T16-59-142	115370	48 ml (115350)용	25ea 6,300
T16-59-145	115372	75 ml (115352)용	20ea 6,400
T16-59-148	115374	100 ml (115354)용	20ea 7,500
T16-59-151	115376	160 ml (115356)용	10ea 3,700

- Non-Autoclavable



Centrifuge Tube, Round Bottom 플라스틱 원심관


- 재질 : 반투명 PP, -10 ~ 135°C 내열, Autoclavable
투명 PC, -135 ~ 135°C 내열, Autoclavable
- Max RCF : 50.000 xg • Tarsons

원심관


Cat. No.	Model	용량	재질	Size (Ø×h)	PK Price ₩
T16-46-168	T542010	16 ml	PP ①	18×100 mm	12ea 12,700
T16-46-176	T542020	50 ml		29×102 mm	12ea 21,100
T16-46-184	T543010	16 ml	PC ②	18×100 mm	12ea 38,800
T16-46-192	T543020	50 ml		29×102 mm	12ea 62,700

원심관 캡, PP Cap ③

Cat. No.	Model	규격	PK Price ₩	CS Price ₩
T16-46-222	T540010	16 ml 원심관용	10ea 3,100	10pk 28,800
T16-46-230	T540020	50 ml 원심관용	10ea 3,800	10pk 36,200



PP



PC



PC Centrifuge Tube, Round Bottom 투명 플라스틱 눈금형 원심관

Cat. No.	Model	용량	Size (Ø×h)	눈금	EA Price ₩
T16-59-191	781029	50 ml	Ø35×99 mm	1 ml	55,700
T16-59-193	781039	100 ml	Ø41×115 mm	2 ml	81,400
T16-59-195	781040		Ø45×98 mm	1 ml	76,400

- 재질 : PC, -135 ~ 135°C 내열, Autoclavable
- Embossed Graduation, 눈금 간격 약 1 mm의 투명한 원심관
- Max. RCF : 20,000 xg • Brand GmbH



Graduated Centrifuge Tube 눈금형 원심관

Cat. No.	Model	규격	PK Price ₩	CS Price ₩
T16-07-060	10-0122	PP Tube	100ea 13,200	10pk 125,400
T16-07-067	10-0123	PE Cap	100ea 4,600	10pk 43,700

- 재질 : PP (원심관), PE (캡) • Autoclavable (Only Tube)
- Max. RCF 3.000 xg • 12 ml 까지 몰드 눈금 표시 (전체 용량 15 ml)
- Size : Ø17.2×h 110.3 mm • Biologix

